

**LUDLUM MODELS 4901P
PANCAKE G-M
HAND AND SHOE MONITOR**

**March 2005
Serial No. 215373 and Succeeding
Serial Numbers**



LUDLUM MEASUREMENTS, INC.

**501 OAK ST., P.O. BOX 810
SWEETWATER, TX 79556
325/235-5494 FAX: 325/235-4672**

STATEMENT OF WARRANTY

Ludlum Measurements, Inc. warrants the products covered in this manual to be free of defects due to workmanship, material, and design for a period of twelve months from the date of delivery. The calibration of a product is warranted to be within its specified accuracy limits at the time of shipment. In the event of instrument failure, notify Ludlum Measurements to determine if repair, recalibration, or replacement is required.

This warranty excludes the replacement of photomultiplier tubes, G-M and proportional tubes, and scintillation crystals which are broken due to excessive physical abuse or used for purposes other than intended.

There are no warranties, express or implied, including without limitation any implied warranty of merchantability or fitness, which extend beyond the description of the face there of. If the product does not perform as warranted herein, purchaser's sole remedy shall be repair or replacement, at the option of Ludlum Measurements. In no event will Ludlum Measurements be liable for damages, lost revenue, lost wages, or any other incidental or consequential damages, arising from the purchase, use, or inability to use product.

RETURN OF GOODS TO MANUFACTURER

If equipment needs to be returned to Ludlum Measurements, Inc. for repair or calibration, please send to the address below. All shipments should include documentation containing return shipping address, customer name, telephone number, description of service requested, and all other necessary information. Your cooperation will expedite the return of your equipment.

**LUDLUM MEASUREMENTS, INC.
ATTN: REPAIR DEPARTMENT
501 OAK STREET
SWEETWATER, TX 79556**

**800-622-0828 325-235-5494
FAX 325-235-4672**

LUDLUM MEASUREMENTS, INC.

501 OAK ST. / P.O. BOX 810

SWEETWATER, TEXAS 79556

Phone: 325-235-5494 800-622-0828 (USA)

Fax: 325-235-4672 E-Mail: ludlum@ludlums.com

Website: <http://www.ludlums.com>**DESIGNER AND MANUFACTURER
OF***Scientific and Industrial
Instruments*

M 4901P Assembly and Disassembly Instructions

The hand detector vertical sections are now shipped detached from the foot detector section. Four (4) screws (8-32X1/2") are used to reattach the vertical sections. A connector is used to distribute power/signals to and from the detectors and main electronics.

✓**Note:** The floor pan is wired such that either upright section may be attached to either side. The floor pan therefore is "non-polarized," and the main electronics will recognize the right and left foot detectors correctly.

Suggested Assembly (Setup) Procedure:

- 1) Carefully unpack the two upright sections and the floor pan section.
- 2) Loosen the four screws located on the end of the foot detector section. Leave the upper two screws in place with about 1/4" of thread showing. Remove the lower two screws.
- 3) Lay one of the uprights (detector face down) on the floor or workbench near the opening on either end of the foot section.
- 4) Look inside the opening for the header that will accept the red plug at the lower end of the upright. Carefully attach the plug to this header.

✓**Note:** The wires should exit the header/plug pointing downward. Make sure the plug is positioned properly (there should be no pins showing on either side of the plug).

- 5) Carefully raise the upright and hang the assembly on the two screws that were left in step 2 above. The upper holes in the ears of the upright are slotted.
- 6) Start the two lower screws and tighten all four of these securely.
- 7) Repeat steps 2 through 6 above for the remaining upright section.
- 8) Attach the power cord and turn the unit ON.
- 9) Check that the unit returns to normal service (**READY LED** will light) after the 60-second update interval has expired.

LUDLUM MEASUREMENTS, INC.

501 OAK ST. / P.O. BOX 810

SWEETWATER, TEXAS 79556

Phone: 325-235-5494 800-622-0828 (USA)

Fax: 325-235-4672 E-Mail: ludlum@ludlums.com

Website: <http://www.ludlums.com>



**DESIGNER AND MANUFACTURER
OF**

*Scientific and Industrial
Instruments*

Suggested Disassembly (Teardown) Procedure (over)

M 4901P Assembly and Disassembly Instructions (continued)

Suggested Disassembly (Teardown) Procedure

- 1) Turn the power OFF to the Model 4901P and remove the power cord from the receptacle.
 - 2) Place the unit on a workbench or other suitable work area.
 - 3) Loosen the four screws holding one of the upright sections.
 - 4) Leave the upper two screws in place and completely remove the lower screws.
 - 5) Carefully lift the upright off and away from the foot section while disconnecting the harness from the floor header.
 - 6) Reinstall the lower two screws and tighten them to prevent loss.
 - 7) Pack the upright sections and the foot detector section well enough to prevent contact with each other and to provide good cushioning.
- ✓**Note:** At least two inches of packing should be provided.

Model 4901P Hand & Shoe Monitor

TABLE OF CONTENTS

1. GENERAL	1
2. SPECIFICATIONS	1
3. SAFETY CONSIDERATIONS AND WARNING MARKINGS	2
4. DESCRIPTION OF CONTROLS AND FUNCTIONS	2
5. ASSEMBLY INSTRUCTIONS	4
5.1 Assembly (Setup) Procedure	4
5.2 Disassembly (Tear-down) Procedure	4
6. SETUP OPERATION	5
6.1 Setup Menu	5
6.1.1 Setup Alarm Menu	5
6.1.2 Setup Background Menu	6
6.1.3 Setup Time Menu	6
6.1.4 Setup Volume Menu	7
6.2 Read Menu	8
6.2.1 Read Alarms Menu	8
6.2.2 Read Time Menu	8
6.2.3 Read Volume Menu	8
6.2.4 Password Menu	9
6.3 Cal Menu	9
6.3.1 Display of Hands Count Data	9
6.3.2 Display of Feet Count Data	9
7. USER OPERATION	10
7.1 Cleaning the Instrument	10
8. COMPATIBLE FIRMWARE VERSIONS	10
9. CALIBRATION PROCEDURE	11
9.1 General	11
9.2 Equipment	11
9.3 Annual Calibration Verification Procedure	11
10. TROUBLE SHOOTING	11
PARTS LIST	13
Amp/HVPS Board, Drawing 436 x 53	13
LED Driver Board, Drawing 420 x 4	14
Ballast Board, Drawing 420 x 155	14
Main Board, Drawing 215 x 60	14
LED Display Board, Drawing 436 x 73	15
Preamplifier Board, Drawing 420 x 47	15
Interconnect Board, Drawing 420 x 178	16
Wiring Diagram, Drawing 420 x 162	16
DRAWINGS AND DIAGRAMS	167

Model 4901P Hand & Shoe Monitor

1. GENERAL

The Model 4901P Beta/Gamma Hand & Shoe Monitor is intended for use as a medium-level beta and gamma contamination monitor. There are four count channels in the standard configuration, monitoring the palms of each hand and the soles of each shoe.

The Model 4901P employs a total of twenty-two pancake Geiger-Mueller (GM)-type detectors, five in each hand detector (palm side only) and six in each foot detector. LED indicators show status and alarm location. The Model 4901P allows parameter updating by viewing the built-in 16-character LCD display. Detector counts, background,

alarm set points and all parameters may be viewed on the LCD display.

Switches at each hand detector initiate an interrogation (both switches must activate). Audible alarm and status change indications are standard.

Features of the Model 4901P include: automatic background accumulate with subtract, password protection of parameters, pushbutton adjustment of the alarm audio volume and simple LED status indicators. All parameters are stored in non-volatile memory, requiring no backup battery.

2. SPECIFICATIONS

- **WEIGHT:** 45 lbs.
- **DIMENSIONS:** 29.5" wide X 15" deep X 40" tall.
- **POWER:** 102-132 VAC, 50/60 Hz, 50 watts maximum.
- **FUSE:** 2 each F-1A, 1 amp, 5 x 20mm, 250 volt.
- **BACKGROUND COUNT RATE:**

HANDS: 200 to 250 cpm.

FEET: 250 to 300 cpm.

- **DETECTOR EFFICIENCY (4 pi):**

HANDS:

12% (4pi) Tc-99

12% (4pi) Cs-137

3% (4pi) C-14

Cs-137, 100 sq.cm. yields 7%.

FEET:

10%(4pi) Tc-99

10% (4pi) Cs-137

3% (4pi) C-14 (✓ **Note: This 1 inch diameter source was placed across feet bars where shielding was minimal.**)

Cs-137, 100 sq. cm. yields 4%.

✓ **Note:** Data taken with 25 to 47 mm disc sources placed directly over pancake tube, except where noted.

- **COUNTING CAPACITY:** 9999 counts per minute.
- **SENSITIVITY:** 85 millivolts nominal.
- **COUNT TIME:** Adjustable from 1 to 99 seconds
- **ALARM HOLD TIME:** Adjustable from 1 to 99 seconds
- **AUDIO:** Piezo speaker with keypad adjustable alarm volume.


Model 4901P Hand & Shoe Monitor

- **BACKGROUND TIME:**
Background accumulate time, adjustable from 1 - 99 seconds.
- **BACKGROUND UPDATE INTERVAL:** Adjustable from 1 to 99 minutes. A background count will take place, if the machine is not in use at every interval specified by this timer.
- **FORCE UPDATE:** Background accumulation will be forced within this interval of time if an automatic accumulation has not been made.
- **ENVIRONMENTAL CONDITIONS:**
indoors or outdoors (clement weather only)
no maximum altitude
temperature range of -20°C to 50°C
maximum relative humidity less than 95% (non-condensing)
mains supply voltage range 85 – 265 VAC
maximum transient voltage of 1500 VAC, Installation Category (Overvoltage Category) II (as defined by IEC 1010-1)
Pollution Degree 1 (as defined by IEC 664)

3. SAFETY CONSIDERATIONS AND WARNING MARKINGS

✓ **NOTE:** The operator or responsible body is cautioned that the protection provided by the equipment may be impaired, if the equipment is used in a manner not specified by Ludlum Measurements, Inc.

As per requirements for CE marking, the Model 4901P may be marked with the following warning symbols, in accordance with EN61010-1:

 **ALTERNATING CURRENT (AC)** (IEC 417, No. 5032) - designates an input receptacle that accommodates a power cord intended for connection to AC voltages. Appears on AC panel.



PROTECTIVE CONDUCTOR TERMINAL (per IEC 417, No. 5019) – designates the central grounding point for the safety ground. Visible inside chassis.



CAUTION (per ISO 3864, No. B.3.1). – designates stability issues. During normal use, the stabilizer legs must be installed to avoid a tip-over of the unit. Without the stabilizer legs, a substantial impact to the front or back of the unit could cause the unit to tip and fall. Appears on AC panel.

☑ **NOTE:** Precautions to be taken during cleaning of the portal monitor are specified in Section 7.1.

4. DESCRIPTION OF CONTROLS AND FUNCTIONS

- **READOUT:** LCD, one line, 16-character alphanumeric display.
- **EXIT Key:** Moves back one menu selection.
- **SPEED TRICK:** Press-and-hold the **EXIT** key to quickly return to the **READY** menu.
- **Increment (Up Arrow “↑”) Button:** Moves up one line in the current menu.

Model 4901P Hand & Shoe Monitor

WITHIN PARAMETER SETUP: A digit increments by one. An on/off parameter toggles to the other state.

- **Decrement (Down Arrow "↓")**

Button: Moves down one line in the current menu.

WITHIN PARAMETER SETUP: A digit decrements by one. An on/off parameter toggles to the other state.

- **SELECT Key:** Selects the current menu choice.

☞ **SPEED TRICK:** If the **SELECT** key is held down while a count channel is being displayed, the alarm level for that channel displays. Releasing the **SELECT** key returns to the count for that channel.

- **SAVE Key:** Recessed pushbutton that saves all parameters to non-volatile memory. This button can only be operated by inserting a small screwdriver or pin through the hole. All of the microprocessor RAM is transferred to flash memory when this button is pushed. Any changes made to variables only change the current microprocessor RAM. If the Model 4901P is turned off prior to saving changes, these changes are lost. To save parameter changes, press the **SAVE** button before turning the Model 4901P off. Upon power-up, the flash memory is loaded into the microprocessor.

- **LEDS:** (Refer to drawings at the back of this manual for LED layout)

- **READY LED:** Must be lit prior to any interrogation.

NOTE: It is possible to begin an interrogation from any setup prompt (when all LED's are on). A count may be started by pressing either of the hand switches, causing the **READY** LED to activate, followed by the **COUNTING** LED.

- **COUNTING LED:** Indicates that a hand count is in progress. Deactivating either of the hand switches prior to expiration of the count time will cause this light to go off and the **SHORT COUNT** LED to come on. When the count is complete, the **CHECK OK** or **ALARM** LED activates. **CHECK OK** lights for 2 seconds or until the hand switches are released. The **SHORT COUNT** LED stays on for the alarm hold time or until a hand switch is reactivated. If no hand switch is sensed during the short count time, then the Model 4901P goes back to the ready state, lighting the **READY** LED.

- **CHECK OK LED:** Indicates that a count has been completed and no alarms were sensed. This LED will stay on until the hand switches are released or for 2 seconds.

- **ALARM:** Indicates that a count has exceeded the alarm set point. The individual **LH, RH, LF or RF** LED lights as soon as alarm is sensed and remain(s) lit for the alarm hold time. The Model 4901P will return to the ready state. The master **ALARM** LED and audible alarm will activate after the end of the count time and the user has removed both hands from the switches. This condition will exist for the duration of the **ALARM HOLD TIME**.

- **SHORT COUNT LED:** Indicates that a count was in progress and the user raised off either hand switch. The **SHORT COUNT** LED stays on for the alarm hold time or until the hand switches are re-activated. A short count resets the count time. If no hand switch is sensed during the short count time, then the Model 4901P goes back to the ready state lighting the **READY** LED.

- **POWER/OK LED:** Indicates that 5Vdc is available on the central processor board.

- **LH, RH, LF, RF LED's:** Indicate which channel in a count has alarmed. These

Model 4901P Hand & Shoe Monitor

light as soon as an alarm is sensed. When the count is complete, the **CHECK OK** will not light. The **ALARM LED** and audio stay on for the alarm hold time, then the Model 4901P goes back to the ready state lighting the

READY LED.

- **Power On/Off:** Switch to turn instrument on and off.

5. ASSEMBLY INSTRUCTIONS

The hand detector vertical sections are shipped detached from the foot detector section. Four (4) screws (8-32X1/2") are used to reattach the vertical sections. A connector is used to distribute power/signals to and from the detectors and main electronics.

Note: The floor pan is wired such that either upright section may be attached to either side. The floor pan therefore is "non-polarized," and the main electronics will recognize the right and left foot detectors correctly.

5.1 Assembly (Setup) Procedure

1. Carefully unpack the two upright sections and the floor pan section.
2. Loosen the four screws located on the end of the foot detector section. Leave the upper two screws in place with about 1/4" of thread showing. Remove the lower two screws.
3. Lay one of the uprights (detector face down) on the floor or workbench near the opening on either end of the foot section.
4. Look inside the opening for the header that will accept the red plug at the lower end of the upright. Carefully attach the plug to this header. Note: The wires should exit the header/plug pointing downward. Make sure the plug is positioned properly (there should be no pins showing on either side of the plug).

5. Carefully raise the upright and hang the assembly on the two screws that were left in step 2 above. The upper holes in the ears of the upright are slotted.

6. Start the two lower screws and tighten all four of these securely.

7. Repeat steps 2 through 6 above for the remaining upright section.

8. Attach the power cord and turn the unit ON.

9. Check that the unit returns to normal service (**READY LED** will light) after the 60-second update interval has expired.

5.2 Disassembly (Teardown) Procedure

1. Turn the power OFF to the M4901P and remove the power cord from the receptacle.
2. Place the unit on a workbench or other suitable work area.
3. Loosen the four screws holding one of the upright sections.
4. Leave the upper two screws in place and completely remove the lower screws.
5. Carefully lift the upright off and away from the foot section while disconnecting the harness from the floor header.
6. Reinstall the lower two screws and tighten them to prevent loss.

7. Pack the upright sections and foot detector section well enough to prevent contact with each other and to provide good

cushioning. **Note:** at least two inches of packing should be provided.

6. SETUP

This section gives instructions on how to use the keys to setup the instrument. Examples of keystroke sequences are given for each parameter. For information on using the instrument to make a radiation check, see section 7.

6.1 Setup Menu

The setup menu has six choices:

- 1- Setup ALARMS MENU
- 2- Setup BACKGROUND MENU
- 3- Setup CAL MENU
- 4- Setup PASSWORD MENU
- 5- Setup TIME MENU
- 6- Setup VOLUME MENU

To change a parameter, access the variable of interest through the setup menus using the **SELECT** and **increment/decrement** “ \uparrow/\downarrow ” keys. Press the **SELECT** key to change the parameter. The cursor becomes visible and blinks on the variable to change. On multiple digit variables, press the **SELECT** key to access the next digit.

☞ **SPEED TRICK:** After changing a parameter, press and hold **SELECT** until a beep is heard. This will quickly exit the setup parameter mode. The setup mode has a blinking cursor.

6.1.1 Set up Alarm Menu

The **SETUP ALARM** menu allows changes to be made to the individual count alarms. All alarm and background values are in units of Counts per Minute.

• INDIVIDUAL ALARMS

The individual channel alarms are lefthand, righthand, leftfoot and rightfoot (LH, RH, LF, and RF).

If the counts are greater than or equal to the count alarm set point for an individual channel during the count time, then the individual alarms LED's (LH, RH, LF, RF) activate. When the count time expires and an alarm is present, the alarm audio sounds and the main **ALARM** LED activates. The alarm will sound for the preset ALARM HOLD TIME.

To access the **SETUP ALARM** menu:

☐ Turn the instrument ON. Wait for **READY** to display on LCD.

☐ Press **SELECT** once, to select the setup menu. **SETUP** menu appears.

☐ Press **SELECT** once to execute the setup menu. **ALARMS** menu appears.

☐ Press **SELECT** once to execute the alarms menu. **LH ALARM XXXX** appears. The **XXXX** is a number between 0 and 9999. This is the current Left Hand Alarm setting.

☐ To change the current setting press **SELECT** to activate the first digit. Use **increment/decrement** “ \uparrow/\downarrow ” to change first digit as desired. Press **SELECT** to activate the second digit. Use **increment/decrement** “ \uparrow/\downarrow ” to change second digit as needed. Press **SELECT** to temporarily save the setting.

☐ From the **LH ALARM XXXX** selection, the increment/decrement keys may be pressed to access further parameters.

Model 4901P Hand & Shoe Monitor

☐ Press the **EXIT** key to exit back to the **ALARMS** menu.

NOTE: Activate the **SAVE** function in order to store all new parameters in non-volatile memory before power down. A small screwdriver, or other object must be used to activate the save feature.

• **LOW BACKGROUND ALARMS**

Set the parameter for **LO BKGND-RH** to a value that would allow detection of a bad detector. For backgrounds near 100 counts in one minute this might be 50. Set the **LO BKGND-LH, LF, and RF** parameter to similar values.

• **HIGH BACKGROUND ALARMS**

Set the high background parameters to preclude nuisance alarms from varying backgrounds. For backgrounds near 100 counts per minute, choose 175. Set both the **LH/RH** and **LF/RF** high background set points.

6.1.2 Setup Background Menu

Access the **SETUP** menu:

- ☐ With **READY** displayed on LCD.
- ☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.
- ☐ Press **SELECT** once again to execute the setup menu. **ALARMS** menu appears.
- ☐ Press **decrement "↓"** once to advance to the **BACKGRND MENU**.
- ☐ Press **SELECT** once to activate menu.
- ☐ Press **SELECT** and use either **increment** or **decrement "↑/↓"** key to toggle the background subtract feature on or off as desired. This will normally be left in the On

position. Activate and exit the on/off prompt by pressing the **SELECT** key one last time.

☐ Press the **decrement "↓"** key to move to the **FORCE UPDATE** interval timer. Press the **SELECT** key to edit this timer as desired. This interval is the maximum time allowed between updates and would normally be set to 15 or 30 minutes. This parameter should be set prior to setting the Update Interval Time or Background Count Time and must always be larger than or equal to either of those (see below). Save and exit this menu item by pressing the **SELECT** key one last time.

☐ Press the **decrement "↓"** key to select the **BKGND UPD INT** timer. This parameter sets the time that will elapse after a hand switch event has ended and a background update takes place. Typical settings are 01 minute. Save and exit this item by pressing the **SELECT** key one time. Note: this parameter must be greater than or equal to the **BKGND TIME** parameter below and less than or equal to the **FORCE UPDATE** parameter above.

☐ Press the **decrement "↓"** key to select the **BKGND TIME**. This is the actual background count time and may be set from 1 to 99 seconds. Longer count times will tend to smooth the background subtract data. Typical count times might be 60 seconds. Note: This number must be less than or equal to the **FORCE UPDATE** and **BKGND UPD INT** parameters as described above.

6.1.3 Setup Time Menu

This menu sets the count time and alarm hold time. The alarm hold time also applies to the **SHORT COUNT LED**.

To access the **SETUP TIME** menu:

- ☐ With **READY** displayed on the LCD.

Model 4901P Hand & Shoe Monitor

☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.

☐ Press **SELECT** once again to execute the setup menu. **ALARMS** menu appears.

☐ Press the **increment** “**↑**” key twice. **TIME MENU** appears.

☐ Press **SELECT** once to execute the setup time menu. **COUNT TIME XX** appears. The XX is a number between 0 and 99 (seconds).

☐ Press **SELECT** to activate the first digit. Use **increment/decrement** “**↑/↓**” to change the first digit. Press **SELECT** to activate the second digit. Use **increment/decrement** “**↑/↓**” to change the second digit. Press **SELECT** to temporarily save parameter.

☐ Use **increment/decrement** “**↑/↓**” to change to the next setting.

☐ Press the **EXIT** key to exit back to the **TIME** menu.

• COUNT TIME

The count time is adjustable between 1 and 99 seconds. This time applies to a count activated by the hand switches. Both of the hand switches must be held down for the duration of the count. If they are not, the **SHORT COUNT LED** activates.

• ALARM HOLD TIME

The alarm hold time is adjustable from 1 to 99 seconds. This time applies to a hand count that has alarmed. If the **ALARM LED** lights, then this light and alarm audio will be held for the alarm hold time. The **SHORT COUNT LED** will also light for this hold time.

NOTE: Remember to press the **SAVE** key in order to store parameters in non-volatile memory prior to power down.

6.1.4 Setup Volume Menu

The volume menu sets only the **ALARM** volume. The Model 4901P emits a beeping sound after various events (mode change, parameter change, etc.). This beeping volume is always at the maximum and is not adjustable.

To access the **SETUP VOLUME** menu:

☐ **READY** is displayed on LCD.

☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.

☐ Press **SELECT** once again to execute the setup menu. **ALARMS** menu appears.

☐ Press the **increment** key once. **VOLUME MENU** appears.

☐ Press **SELECT** once to execute the setup volume menu. **ALARM VOLUME XXX** appears. The XXX is a number between 0 and 255. This variable sets from 255 (lowest level) to 000 (maximum level). Any audio alarm uses this volume set point. The beep audio is not affected by this setting.

☐ Press **SELECT** to activate the first digit. Use **increment/decrement** “**↑/↓**” to change the first digit. Press **SELECT** to activate the second digit. Use **increment/decrement** “**↑/↓**” to change the second digit. Repeat for third digit. Press **SELECT** to save.

☐ Press the **EXIT** key to exit back to the **VOLUME** menu.

NOTE: Remember to press the **SAVE** key in order to store any changed parameters in non-volatile memory prior to power down.

Model 4901P Hand & Shoe Monitor

6.2 Read Menu

The read menu has three choices:

- 1- Read Alarms Menu
- 2- Read Time Menu
- 3- Read Volume Menu

The read menu accesses the same menu structure as the Setup Menu. **However, no variables may be changed from the read menu.**

6.2.1 Read Alarms Menu

To access the READ ALARMS menu:

- ☐ Turn the instrument ON. Wait for **READY** to display on LCD.
- ☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.
- ☐ Press increment key “**↑**” once. **READ** menu appears.
- ☐ Press **SELECT** once to execute the read menu. **ALARM** menu appears.
- ☐ Press **SELECT** once to execute the alarms menu. **GLOBAL ALARM XX** appears. The XX is a number between 0 and 99.
- ☐ Use the increment/decrement “**↑/↓**” keys to change to the next alarm channel.
- ☐ Press the **EXIT** key to exit back to the **ALARMS** menu.

6.2.2 Read Time Menu

This menu reads all of the time parameters of the Model 4901P. The user cannot change these values from this menu.

To access the READ TIME menu:

- ☐ Turn instrument ON. Wait for **READY** to display on LCD.
- ☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.
- ☐ Press **decrement** key “**↓**” once. **READ** menu appears.
- ☐ Press **SELECT** once to execute the read menu. **ALARMS** menu appears.
- ☐ Press decrement key “**↓**” once. **TIME** menu appears.
- ☐ Press **SELECT** once to execute the time menu. **COUNT TIME XX** appears. The XX is a number between 0 and 99.
- ☐ Use the increment/decrement “**↑/↓**” keys to change to other time parameters.
- ☐ Press the **EXIT** key to exit back to the **TIME** menu.

6.2.3 Read Volume Menu

This menu reads all of the volume parameters of the Model 4901P. The user cannot change these values from this menu.

To access the READ VOLUME menu:

- ☐ Turn the instrument ON. Wait for **READY** to display on LCD.
- ☐ Press **SELECT** once to select the setup menu. **SETUP** menu appears.
- ☐ Press decrement key “**↓**” once. **READ** menu appears.
- ☐ Press **SELECT** once to execute the read menu. **ALARMS** menu appears.

Model 4901P Hand & Shoe Monitor

☐ Press decrement “↓” key twice. VOLUME menu appears.

☐ Press **SELECT** once to execute the time menu. ALARM VOLUME XXX appears. The XXX is a number between 0 and 255.

☐ Use the increment/decrement “↑/↓” keys to change to other parameters.

☐ Press the **EXIT** key to exit back to the VOLUME menu.

6.2.4 Password Menu

This menu sets the password and whether the password is On or Off.

To access the PASSWORD menu:

☐ With READY displayed on the LCD.

☐ Press **SELECT** once to select the setup menu. SETUP menu appears.

☐ Press **SELECT** once to execute the setup menu. ALARMS menu appears.

☐ Press the increment or decrement “↑/↓” keys until the PASSWORD menu appears.

☐ Press **SELECT** once to execute the password on/off menu. PASSWORD: XXX appears. The XXX is either ON or OFF.

☐ Press **SELECT** to change the password status. Use increment/decrement “↑/↓” to change to either ON or OFF. Press **SELECT** to temporarily save parameter.

☐ Use increment/decrement “↑/↓” to change to the next setting. ENTER PASS: XXXX appears.

☐ To reset the PASSWORD to 0000, hold down the SAVE key while turning on the instrument.

☐ Press **SELECT** to activate the first digit. Use increment/decrement “↑/↓” to change the first digit. Press **SELECT** to activate the second digit. Use increment/decrement “↑/↓” to change the second digit. Repeat for third and fourth digit. Press **SELECT** to save.

☐ Press the **EXIT** key to exit back to the TIME menu.

NOTE: Press the SAVE key in order to store parameters in non-volatile memory prior to power down.

6.3 Cal Menu

The Cal menu has two choices:

6.3.1 Display of Hands Count Data

Selecting this mode provides a one second updating display of the current count from the hand detectors (in counts per second). This mode is used for setting or checking the threshold level and as a general diagnostic using a pulser or source counts from the detectors.

6.3.2 Display of Feet Count Data

Selecting this mode provides fast, one second updating display of the current count from the feet detectors (in counts per second). This mode is used for setting or checking the threshold level and as a general diagnostic using a pulser or source counts from the detectors.

7. USER OPERATION

This section gives instructions on how to use the instrument to make a radiation check. For information on Parameter Setup, see Section 5.

A count starts when both of the hand switches are held down. If the LCD was in a SETUP menu, then the LCD returns to the READY menu and a normal count will take place. If the LCD was in READ COUNTS menu the LCD will remain in this menu and the interrogation will proceed normally. Note: When monitoring counts via the CAL MODE, an interrogation will not be available (the READY LED will be extinguished).

Prior to operation, the monitor must be allowed to update the background count. This mandatory update occurs just after power-up and after expiration of the Force Update interval timer. New background count data is compared to the low and high background set points that have been programmed into the unit. If the set points have been exceeded, an alarm is given (check individual LED's for offending channel) and the unit returns to updating background.

In order to make a radiation check, follow the steps below.

☐ The green **READY** light must be lit in order to use the instrument.

☐ Step up and position both hands over the detectors.

☐ Place palms flat against the bottom screen and push inward until the green **COUNTING** light turns on.

☐ The yellow **SHORT COUNT** light will turn on if the hands are removed before the count is complete.

☐ Once count is complete, the green **CHECK OK** light or the red **ALARM** light will turn on. Smaller red lights will turn on with the **ALARM** light to indicate the location of the alarm.

☐ Remove hands and step off instrument.

7.1 Cleaning the Instrument

The Model 4901P may be cleaned with a damp cloth (using only water as the wetting agent). Do not immerse the instrument in any liquid. Observe the following precaution when cleaning:

Turn instrument OFF and disconnect the instrument power cord.

8. COMPATIBLE FIRMWARE VERSIONS

FIRMWARE- A computer program loaded into permanent memory of the instrument. This hardware (memory) cannot be changed in its user environment.

This manual works with instrument firmware versions:

M4901P: 420-03-N01

The firmware number displays when the instrument is first turned on or may be viewed through the diagnostic menu.

9. CALIBRATION PROCEDURE

9.1 General

The Model 4901P was set up for 80 mV sensitivity and 900 Vdc operation for G.M. type detectors.

9.2 Equipment

1. Ludlum Model 500 Pulser or equal
2. High Impedance voltmeter for high voltage measurements (10 megohm)
3. 8 to 15 volt DC power supply with modular connection (pin 2 is positive and pin 3 is ground) polarity protected

9.3 Annual Calibration Verification Procedure

Calibration of the Model 4901P is accomplished by checking the threshold level

at each preamplifier board (LMI #5436-040) located on each detector.

The design threshold level is 80 mV and operating high voltage is approximately 900 Vdc.

☐ Using a clip lead cable, connect the Model 500 Pulser to the detector ballast board and apply power to the board.

☐ Sweep pulser amplitude for a negative leading edge 70 to 90 mV pulse and confirm counter turn on at 80 mV +/- 5 mV. If necessary adjust R1 (THS) until pulses just appear.

☐ Check for 900 Vdc +/- 10 V at the detector ballast board. If necessary, adjust R4 (HV ADJ) for 900 Vdc at the ballast board input.

10. TROUBLE SHOOTING

The block diagram of the M4901P can be thought of as four detectors connected to a multi-counter MAIN board. All detectors operate from a single, high voltage power supply (HVPS). This supply is located on the Main Electronics chassis just below the Main board. The count data appears at this Main board as 5-volt digital pulses. These pulses are generated on the preamplifier board at each of the four detectors. Calibration is performed on each detector and consists of setting the lower level threshold or discriminator (LLD) and setting the HV bias to the proper operating point.

The User LED board presents status information to the user via a serial data stream from the Main controller board. This serial data is placed into two drivers that directly drive the LED's.

The Main control board also sends data to the LCD display. The LCD is intended for setup purposes as well as diagnostics. It is not necessary for the user to view the LCD screen under normal conditions. Count data can be reviewed in the display if desired.

The "pancake" G.M. detectors used in this model are simple in application but can cause headaches when "ganged" in parallel as in the M4901P. One bad detector can cause the entire unit to become noisy, due mainly to the use of the single HVPS. Normally, only one detector becomes noisy and the culprit can be found in that particular array. A quick visual check may reveal the bad detector. Inspect the thin membrane cover of each of the pancakes to see if one of them has lost its gas. The membrane will look loose or wrinkled and when touched (carefully) will make a

Model 4901P Hand & Shoe Monitor

crackling sound. This one will definitely need replacement. If you find no broken membranes and you are in a relatively quiet area, you can listen to each tube for the one that is noisy. Each event in the tube is an avalanche of charge (a spark) so they can be heard rather easily, provided you have adequate HV bias. A single bad probe can pull the HV bias down and prevent all others from working.

As a last resort the detector array in question will have to be removed and each

detector signal wire unplugged until the offending pancake has been located. The signal wires have a connector on one end to facilitate fast, no-solder removal.

There are no batteries required for parameter storage during power down. All parameters are saved in Flash memory when the Store button is pressed. Press store anytime you change parameters and wish them to be used from then on. If you do not press store, the old values will reappear after the next power down and up cycle.

CAUTION:
TO AVOID ELECTRICAL SHOCK, ENSURE
THAT THE INSTRUMENT IS OFF FOR AT LEAST
ONE MINUTE BEFORE TOUCHING THE CONNECTIONS.

Model 4901P Hand & Shoe Monitor

PARTS LIST

Ref. No.	Description	Part No.	Ref. No.	Description	Part No.	
Model 4901P Hand & Shoe			Q131	MTD10N05E	05-5839	
<hr/>			Q211	PQ20VZ51	05-5863	
UNIT	Completely Assembled M4901P	48-3009	Q221	MMBT3904T	05-5841	
	Hand & Shoe Monitor		Q222	MMBT4403LT1	05-5842	
<hr/>			Q223	TIP120	05-5782	
HVPS Board, Drawing 436 x 53			<hr/>			
BOARD	Assembled HVPS	5436-042	CONNECTORS			
<hr/>			J130	CONN		
CAPACITORS			P3,P5,P6	RAPC712 93F7715	13-8445	
C001	10μF 20V	04-5655		CONN-640456-2		
C002	1μF 35V	04-5656	P4	MTA100	13-8073	
C011-C014	0.0047μF 3KV C	04-5547		CONN-640456-4		
C021-C023	0.0047μF 3KV C	04-5547		MTA100	13-8088	
C024	0.0027μF 3KV C NPO	04-5520	RESISTORS			
C031	0.0027μF 3KV C NPO	04-5520	R001	2.21K 1/8W 1%	12-7835	
C101	1μF 35V	04-5656	R002	3.32K 1/8W 1%	12-7870	
C102	10μF 20V	04-5655	R003	2.21K 1/8W 1%	12-7835	
C111	0.0047μF 3KV C	04-5547	R011	475K 1/8W 1%	12-7859	
C112	0.01μF 50V X7R	04-5664	R012	1 GIG-OHM FHV-1 2%	12-7686	
C113	0.1μF 50V X7R	04-5663	R013	TRMR-1 MEG	09-6911	
C114	0.01μF 50V X7R	04-5664	R014-R015	100 K 1/4W 5%	10-7023	
C121	100pF 3KV 30GAT10	04-5532	R021	100 K 1/4W 5%	10-7023	
C122	0.0047μF 3KV C	04-5547	R111	1M 1/8W 1%	12-7844	
C123	100pF 100V COG	04-5661	R112-R113	1 GIG-OHM FHV-1 2%	12-7686	
C124	0.1μF 50V X7R	04-5663	R114	10 MEG 1/4W 5%	12-7955	
C128	0.1 μF 16V	04-5701	R115	1M 1/8W 1%	12-7844	
C131	68μF 6.3V	04-5654	R116	TRMR-1 MEG	09-6911	
C211	47μF 10V	04-5666	R117	1K 1/8W 1%	12-7832	
C212	0.0022μF 50V COG	04-5676	R121	1M 1/8W 1%	12-7844	
C213	47μF 10V	04-5666	R123	432K 1/8W 1%	12-7874	
C214	10μF 20V	04-5655	R124	33.2K 1/8W 1%	12-7842	
C221	10μF 20V	04-5655	R125	182K 1/8W 1%	12-7860	
C231	0.1μF 50V X7R	04-5663	R126	1K 1/8W 1%	12-7832	
C311	1μF 35V	04-5656	R127	4.75K 1/8W 1%	12-7858	
<hr/>			R201	7.5K 1/8W 1%	12-7847	
DIODES			R211	100K 1/8W 1%	12-7834	
CR021-CR022	1N4007	07-6274	R212	165K 1/8W 1%	12-7877	
CR031-CR032	1N4007	07-6274	R213	22.1K 1/8W 1%	12-7843	
CR101	1N5817	07-6290	R214	1.27K 1/8W 1%	12-7902	
DS001	LED-HLMP 3502	07-6280	R215	33.2K 1/8W 1%	12-7842	
DS002	LED-HLMP 3000	07-6288	R221	22.1K 1/8W 1%	12-7843	
<hr/>			R222	4.75K 1/8W 1%	12-7858	
TRANSISTORS			R223	1K 1/8W 1%	12-7832	
Q001	2N7002L	05-5840	R224	TRMR-10K 3269W1-103	09-6931	
Q002	PQ05SZ11 5V	05-5858	R225	18.2 K 1/8W 1%	12-7968	
Q121	2N7002L	05-5840	R226	10K 1/8W 1%	12-7839	
			R227	1K 1/8W 1%	12-7832	
			R228	10K 1/8W 1%	12-7839	
			R330	TRMR-10K 64W103	09-6787	

Model 4901P Hand & Shoe Monitor

Ref. No.	Description	Part No.	Ref. No.	Description	Part No.	
•	INTEGRATED CIRCUITS		•	RESISTORS		
U111	TLC27M7ID	06-6292	R1-R2	100 K 1/4W 5%	10-7023	
U112	TLC372ID	06-6290	R3-R8	3.3 MEG 1/4W 5%	10-7044	
U121	ICM7555CBA	06-6300	•	CONNECTORS		
U211	LT1054CS	06-6315		P1	CONN-640456-3	13-8081
U221	LM285M-1.2	05-5845			MTA100	
U222	LMC7111BIM5	06-6410				
•	TRANSFORMERS		Main Board, Drawing 215 x 60			
T121	XFMR-M 416-3 HV	4275-145	BOARD	Assembled Main	5215-087	
•	MISCELLANEOUS		•	CAPACITORS		
10 EA.	CLOVERLEAF	18-8771	C101	68μF 6.3V	04-5654	
	RECPT-01106809-000		C201	68μF 6.3V	04-5654	
<u>LED Driver Board, Drawing 420 x 4</u>			C211	0.1μF 50V X7R	04-5663	
BOARD	Assembled LED Driver	5420-005	C231	0.01μF 50V X7R	04-5664	
			C301	2700μF 35V E	04-5621	
•	INTEGRATED CIRCUITS		C311	27pF 100V COG	04-5658	
U140-U141	SN75512	06-6369	C312	27pF 100V COG	04-5658	
•	RESISTORS		C501	68μF 6.3V	04-5654	
R148	200 OHM	10-7006	C502	0.1μF 50V X7R	04-5663	
R149-R151	10k	10-7016	C503	10μF 20V	04-5655	
•	RESISTOR NETWORKS		C504-506	0.1μF 50V X7R	04-5663	
RN142-RN144	150 OHM	12-7741	C601	10μF 20V	04-5655	
•	CONNECTORS		C602	4.7μF 20V	04-5653	
P23	CONN-640456-5	13-8057	C603	10μF 20V	04-5655	
	MTA100		C611	4.7μF 20V	04-5653	
<u>BALLAST BOARD, Drawing 420 X 155</u>			C701	0.1μF 50V X7R	04-5663	
BOARD	Assembled Ballast	5420-158	C711	0.1μF 50V X7R	04-5663	
			•	DIODES		
•	CAPACITORS		CR101-103	CXSH-4 EB33	07-6358	
C1	0.0047μF 3KV C	04-5547	•	TRANSISTORS		
			Q211	MMBT4403LT1	05-5842	
			Q401	2N7002L	05-5840	
			Q402	MMBT4403LT1	05-5842	
			Q501	MMBT3904T	05-5841	
			•	CONNECTORS		
			P14	CONN-640456-2	13-8073	
				MTA100		
			P15	CONN-640456-6	13-8095	
				MTA100		
			P16	CONN-640456-3	13-8081	
				MTA100		
			P17	CONN-640456-5	13-8057	
				MTA100		

Model 4901P Hand & Shoe Monitor

Ref. No.	Description	Part No.	Ref. No.	Description	Part No.
P18	CONN-640456-3		S121	92-851.342	08-6726
	MTA100	13-8081	S211	92-851.342	08-6726
P19	CONN-1-640456-1		S221	92-851.342	08-6726
	MTA100	13-8059	S321	92-851.342	08-6726
P20	CONN-1-640456-4		5 EA.	92-960-0 MNT FLANGE	08-6727
	MTA100	13-8141			
•	RESISTORS		•	VOLTAGE REGULATORS	
R031	4.75K 1/8W 1%	12-7858	VR201	LT1129CQ-5	06-6372
R111	100K 1/8W 1%	12-7834	•	RESISTOR NETWORKS	
R131	2.21K 1/8W 1%	12-7835			
R1310	100K 1/8W 1%	12-7834	RN031	NETWORK-4.7 K	12-7918
R132-R139	2.21K 1/8W 1%	12-7835	RN121	NETWORK-4.7K 8P SIP	12-7706
R211-R212	10K 1/8W 1%	12-7839	RN331	NETWORK-4.7 K	12-7918
R231	100K 1/8W 1%	12-7834	RN421	NETWORK-22 K	12-7917
R331	22.1K 1/8W 1%	12-7843	•	CRYSTALS	
R401	10K 1/8W 1%	12-7839			
R402	10 OHM 1/8W 1%	12-7836	Y311	MICRO 6.144 MHZ	01-5262
R403	10K 1/8W 1%	12-7839	•	TRANSFORMERS	
R431	10K 1/8W 1%	12-7839			
R501	10K 1/8W 1%	12-7839			
R502	10 MEG 1/4W 5%	12-7955	T401	M 177 AUDIO	4275-083
R503	73.2K 1/8W 1%	12-7895	•	MISCELLANEOUS	
R504	10K 1/8W 1%	12-7839			
R505	82.5K 1/8W 1%	12-7849			
R506	1M 1/8W 1%	12-7844	*	SOCKET-822276-1 44P	06-6293
R507	8.25K 1/8W 1%	12-7838			
R508	10K 1/8W 1%	12-7839			
R701	TRMR-5K 3269W1-502	09-6918			
•	INTEGRATED CIRCUITS			LED Display Board, Drawing 420 x 73	
			BOARD	Assembled LED Display	5420-097
U121	NETWORK-22K		•	LEDS	
	DIP 14 PIN	12-7577			
U122	LTC1045CN	06-6371	CR110-CR112	LED-E121 GREEN	07-6310
U131	TLC372ID	06-6290	CR113	LED-E176 RED JUMBO	07-6362
U211	X24C02S8I	06-6299	CR114	LED-E120 YELLOW	07-6309
U231-U233	TLC372ID	06-6290	CR115	LED-E121 GREEN	07-6310
U311	N87C51FC	06-6303	CR125-CR130	LED-E112 RED	07-6390
U331	TLC372ID	06-6290	•	MISCELLANEOUS	
U421	CD74HC573M	06-6298			
U431	N82C54	06-6309			
U501	LM358D	06-6312	P12	CONN-CJ50-36B-10	13-8370
U502	LM285M-2.5	06-6291			
U511	CXK581000AM-70LL	06-6385		Preamplifier Board, Drawing 436 x 47	
U531	N82C54	06-6309			
U601	MAX232CSE	06-6382	BOARD	Assembled Preamplifier	5436-040
U611	CD74HC138M	06-6339	•	Capacitors	
U612	CD74HC00M	06-6308			
U631	N82C54	06-6309	C001	0.01μF 50V X7R	04-5664
U711	CD74HC08M	06-6313	C002	10μF 20V	04-5655
•	SWITCHES		C101	10μF 20V	04-5655
			C103	0.001μF 100V COG	04-5659
S111	92-851.342	08-6726			

Model 4901P Hand & Shoe Monitor

Ref. No.	Description	Part No.	Ref. No.	Description	Part No.
C104	0.01μF 50V X7R	04-5664	Interconnect Board, Drawing 420 x 178		
C105	10μF 20V	04-5655			
C106	1.0μF 16V C	04-5701	BOARD	Assembled Interconnect	5420-178
C108	1μF 35V	04-5656			
C109	10μF 20V	04-5655	•	CONNECTORS	
C201	10μF 20V	04-5655	P1	CONN-1-640457-0	
C203	100pF 3KV 30GAT10	04-5532		MTA100-RA	13-8168
•	TRANSISTORS		P2	CONN-1-640456-0	
Q101	MMBT3904T	05-5841		MTA100	13-8066
•	RESISTORS		Wiring Diagram, Drawing 420 x 162		
R001	4.75K 1/8W 1%	12-7858	•	SWITCHES	
R002	100K 1/8W 1%	12-7834	S1	DM62J12S205PQ	08-6715
R003	100 OHM 1/8W 1%	12-7840	S2-S3	BZ-2RD-A2-MICRO	08-6538
R004	100K 1/8W 1%	12-7834	•	TRANSFORMER	
R005	1K 1/8W 1%	12-7832	T1	XFMR-CFP302 115/230V	
R101-R102	47.5 OHM 1/8W 1%	12-7966			22-9908
R104	5.62K 1/8W 1%	12-7871	•	CONNECTORS	
R105	4.75K 1/8W 1%	12-7858	J1	CONN-640456-2	
R106	1.27K 1/8W 1%	12-7902		MTA100	13-8073
R107-R108	2.37K 1/8W 1%	12-7861	J2	CONN-640456-4	
R109	1K 1/8W 1%	12-7832		MTA100	13-8088
R201	1K 1/8W 1%	12-7832	J4	CONN-1-640456-4	
R202-R203	47.5K 1/8W 1%	12-7872		MTA100	13-8141
R304	TRMR-10K 64W103	09-6787	J14	CONN-640456-2	
R1010	22.1K 1/8W 1%	12-7843		MTA100	13-8073
R1011	100 OHM 1/8W 1%	12-7840	J17	CONN-640456-5	
•	INTEGRATED CIRCUITS			MTA100	13-8057
U001	TLC372ID	06-6290	J18	CONN-640456-3	
U101	CA3096M	06-6288		MTA100	13-8081
•	CONNECTORS		J19	CONN-1-640456-1	
P1	CONN-640456-2			MTA100	13-8059
	MTA100	13-8073	J23	CONN-640456-5	
P2	CONN-640456-4			MTA100	13-8057
	MTA100	13-8088	•	MISCELLANEOUS	
•	INDUCTOR		DSO1	UNIMORPH	
L101	INDUCTOR-TKS1245	21-9699		TEC-3526-PU	21-9251

DRAWINGS AND DIAGRAMS

PANEL AND PLATE ASSEMBLY DRAWINGS

Main Chassis Front Panel, Drawing No. 420 x 171

Front Panel LED, Drawing No. 420 x 170

SCHEMATICS AND COMPONENT LAYOUTS

HVPS Board, Drawing No. 436 x 53

HVPS Board Component Layout, Drawing No. 436 x 54

LED Display Driver Board, Drawing No. 420 x 4

LED Display Driver Board Component Layout, Drawing No. 420 x 89

Detector Ballast Board, Drawing No. 420 x 155

Detector Ballast Board Component Layout, Drawing No. 420 x 156

Main Board, Drawing No. (2 sheets) 215 x 60

Main Board Component Layout, Drawing No. 215 x 103

LED Display Board, Drawing No. 420 x 73

LED Display Board Component Layout, Drawing No. 420 x 92

Preamplifier Board, Drawing No. 436 x 47


Preamplifier Board Component Layout, Drawing No. (2 sheets) 436 x 48

Interconnect Board, Drawing No. 420 x 178

Interconnect Board Component Layout, Drawing No. 420 x 179


Wiring Diagram, Drawing 420 x 162

REV #	ALTERATIONS	DATE	BY
	VALID	07-22-98	TJR



LUDLUM
MEASUREMENTS, INC.
SWEETWATER, TEXAS

MODEL 4901P
HAND & SHOE MONITOR



•

• •

•


SAVE

↑

↓


XXXXXX


DATE	REV	DATE	REV
07-22-98	07-22-98	07-22-98	07-22-98
TJR	TJR	TJR	TJR
TITLE: M 4901P MAIN ELEC. PANEL		SHEET 420	
LUDLUM MEASUREMENTS, INC.		SHEET 171	



REV #	ALTERATIONS	DATE	BY
	VALID	07-22-98	TJR


DATE	07-22-98	07-22-98	07-22-98
TJR	07-22-98	07-22-98	07-22-98
TITLE	M 4901P LED DISPLAY		
LUDLUM	MEMBER, INC.	420	170
IN THE FIELD	420	170	170
REVISIONS	420	170	170





LUXUM REPAIRMENTS INC.	
DATE	06/09/84
CHK	06/09/84
DESIGN	06/09/84
APPD	06/09/84
NEXT	HIGHER ASSY.
SIZE	22-58
MODEL	4901/4901-1/52
SERIES	490
SHEET	4
DATE	07/12/87
REV	02-AU-08
REV	02-AU-08

LUXUM REPAIRMENTS INC.	
TITLE: LED DISPLAY DRIVER	
BOARD	5420-065
SIZE	22-58
MODEL	4901/4901-1/52
SERIES	490
SHEET	4
DATE	07/12/87
REV	02-AU-08
REV	02-AU-08





LUDLUM MEASUREMENTS INC. SHEETWATER, TX.	
DR	RDS07JUN94
CHK	CAR 22 02 94
DISC	RDS07JUN94
APP	BS 2-22-94
07:44:11	22-JUL-98
COMP PASTE	COMP MASK
ISDR PASTE	ISDR MASK
TITLE: LED DRIVER	
BOARD#	5420-005
MODEL	BS420005
SERIES	4901/52
SHEET	89
OUTLINE	W




P.O. Box 810
501 Oak Street
Sweetwater, Texas 79556
U.S.A. 1-800-622-0628

Drawn: CKB	17-JUL-00	Title: SIX DETECTOR BALLAST
Design: ROS	17-JUL-00	Model: 4901P
Check: <i>ROS</i>	11-17-02	Board: 5420-159
Approved: <i>ROS</i>	10-27-00	Sheet: 1 of 1
08/34/25	27-Oct-2000	Rev: 2.0
T:\PROJECTS\1000\1000PA0007.Dwg		Document: 2015 Rev: 24/01/15 SGP
Series		420 155







Drawn: CKB	17-JUL-00	Title:			
Design: RDS	17-JUL-00	SIX DETECTOR BALLAST			
Check: P.W.	10-27-00	Model: 4901P			
Approve: RSS	10-27-00	Board#: 5420-158			
Layer:	Top Overlay		Rev. 1.0	Series	Sheet
Mech.1	MD:		SCALE: 1.00	420	156
Mech.2					
Mech.3	08:34:36				
Mech.4	27-Oct-2000				




LUDLUM MEASUREMENTS, INC.		TITLE: MAIN BOARD	
DATE	81-MAY-86	SIZE	MODEL
DR OF	CKR	BOARD	5215-087
REV	1	SIZE	C
APPD	PSS	SERIES	215
BY	87-JAN-86	SHEET	68
BY	87-JAN-86	SHEET	2




UPDATED	-	LUXURY MEASUREMENTS INC.			
BR CIG	08/19/95	TITLE: MAIN BOARD			
CHK FUL	7-22-98				
ISSN RGS	08/16/95	BOARD 5215-987			
APPR	7-22-98	SIZE	MODEL	SERIES	SHEET
NEXT HIGHER ASSY,				215	60
RT:38:12	22-AUG-98	S2	S871587	SHEET 2	OF 2




LUDLUM MEASUREMENTS INC.		SHEETWATER, TX.	
DR	RDS	14SEP94	TITLE: LED DISPLAY BOARD
CHK	CR	23-JUL-96	BOARD# 5420-897
			BS420897
DSON	RDS	14SEP97	MODEL 4901-1
			SERIES 420
APP	BS	22-98	SHEET 92
97:52:149	22-JUL-98	COMP ARTWORK	SLDR ARTWORK
		COMP PASTE	COMP TASK
		SLDR PASTE	SLDR TASK





UPDATED	-	LUDLOW MEASUREMENTS INC.			
DR RDS	07/07/94				
CHK RDS	22-300-13	TITLE: M4981-1 LED DISPLAY			
USN RDS	SEP94	BOARD: 5438-B97			
APPD	22-300-13	SIZE: MODEL 4981-1			
NEXT HIGHER ASSY:		SERIES 428		SHEET 73	
07/17/97	22-300-13	SHEET 73		OF 73	



LUDLUM MEASUREMENTS INC.			
REV	28-AUG-97	TITLE	PREAMP
CHK	7-22-98	BOARD	5408-040
DES	17-JUL-97	SIZE	C
APP	255	MODEL	4901P
NEXT HIGHER ASSY.		SHEET	47
22-AUG-98		OF	47




 LUDLUM MEASUREMENTS INC. SHEETWATER, TX.	
DR	ACF 17-JUL-97 TITLE: PREAMP
CHK	CRB 22-JUL-97 BOARD: 5438-040 BS438040
DESN RDS	17-JUL-97 MODEL 4901P SERIES 438 SHEET 48
APP	BS 2-2-78 COMP. ARTWORK <input type="checkbox"/> SLDR. ARTWORK <input type="checkbox"/>
01:01:97	17-JUL-97 COMP. OUTLINE <input type="checkbox"/> SLDR. OUTLINE <input type="checkbox"/>
COMP. PASTE	<input type="checkbox"/> COMP. MASK <input type="checkbox"/> SLDR. PASTE <input type="checkbox"/> SLDR. TASK <input type="checkbox"/>



LUDLUM MEASUREMENTS, INC.
 PO Box 810
 501 Oak Street
 Sweetwater, Texas 75556
 U.S.A. 1-800-622-0828

Drawn: MO 05-JAN-2000 Title: INTERCONNECT BOARD
 Check: *DM* 10-27-00 Board: 5420-178
 Approv: *DM* 10-27-00 Rev: 1 of 1
 08-38-15 27-Oct-2000 Rev: 1 of 1
 X:\Photos\JAW001\JAW001.Dwg (Source: JAW001.Dwg) 10/27/2000 1:18 PM

Sheet 420 of 178



Drawn: MG		06-JAN-2000	Title:		
Design: RDS		06-JAN-2000	INTERCONNECT BOARD		
Check: PW.		10-27-00	Model: M4901P		
Approve: BS		10-27-00	Board#: 5420-178		
Layer:	Top Overlay		Rev: 1.0	Series	Sheet
Mech.1	MD:				
Mech.2	08:42:05		SCALE: 1.00	420	179
Mech.3	27-Oct-2000				
Mech.4					
bs420178.pcb			bs420178.pcb		

