

**LUDLUM MODEL 2350-1
DATA LOGGER**

Calibration Routines, Parts List, and Schematics

April 2022

**Serial No. 331990 and Succeeding
Serial Numbers**

LUDLUM MODEL 2350-1 DATA LOGGER

Calibration Routines, Parts List, and Schematics

April 2022

**Serial No. 331990 and Succeeding
Serial Numbers**



LUDLUM MEASUREMENTS, INC
501 OAK STREET, P.O. BOX 810
SWEETWATER, TEXAS 79556
325-235-5494, FAX: 325-235-4672

STATEMENT OF WARRANTY

Ludlum Measurements, Inc. warrants the products covered in this manual to be free of defects due to workmanship, material, and design for a period of twelve months from the date of delivery. The calibration of a product is warranted to be within its specified accuracy limits at the time of shipment. In the event of instrument failure, notify Ludlum Measurements to determine if repair, recalibration, or replacement is required.

This warranty excludes the replacement of photomultiplier tubes, G-M and proportional tubes, and scintillation crystals which are broken due to excessive physical abuse or used for purposes other than intended.

There are no warranties, express or implied, including without limitation any implied warranty of merchantability or fitness, which extend beyond the description of the face there of. If the product does not perform as warranted herein, purchaser's sole remedy shall be repair or replacement, at the option of Ludlum Measurements. In no event will Ludlum Measurements be liable for damages, lost revenue, lost wages, or any other incidental or consequential damages, arising from the purchase, use, or inability to use product.

RETURN OF GOODS TO MANUFACTURER

If equipment needs to be returned to Ludlum Measurements, Inc. for repair or calibration, please send to the address below. All shipments should include documentation containing return shipping address, customer name, telephone number, description of service requested, and all other necessary information. Your cooperation will expedite the return of your equipment.

**LUDLUM MEASUREMENTS, INC.
ATTN: REPAIR DEPARTMENT
501 OAK STREET
SWEETWATER, TX 79556**

**800-622-0828 325-235-5494
FAX 325-235-4672**

TABLE OF CONTENTS

1. INITIAL INSTRUMENT CALIBRATION	4
1.1 Equipment Required	4
1.2 Instrument Preparation	4
1.3 Threshold and Window Calibration	4
1.4 High Voltage Calibration	5
1.5 Overload Calibration	5
1.6 Alarm Setpoints	6
2. CALIBRATION ROUTINES	7
Source Parameters and Typical Values for Dead Time and Calibration Constants .	7
2.1 High Voltage Ramp Routine	8
2.2 Calibration Constant and Dead Time Calibration with Background Subtract	10
2.3 Single Point Calibration Constant Routine with Background Subtract	12
2.4 Two Source Dead Time Calibration Routine with Background Subtract	14
3. THEORY OF OPERATION	16
3.1 Amplifier/Power Supply Board	16
3.2 Central Processor Board	18
3.3 I/O Interface Board	19
3.4 64k Memory Expander Board	20
4. PARTS LIST, COMPONENT LAYOUT, AND SCHEMATICS	21
Amplifier / Power Supply Board	21
Central Processor Board	24
64k Memory Expander Board	26
I /O Interface Board.....	27
Backplane Board.....	28
List of Drawings.....	28

MODEL 2350-1 Data Logger

1. INITIAL INSTRUMENT CALIBRATION

1.1 Equipment Required

All instruments used in calibrating the Model 2350-1 must be calibrated by standards traceable to the National Institute of Standards and Technology and must have a current calibration label attached.

1. Ludlum Model 500 Pulser.
2. High Impedance voltmeter with at least 1000 megohms meter resistance and an accuracy of 0.5%.
3. Overload simulator (1000 megohm resistor).



1.2 Instrument Calibration

The following procedures will calibrate the electronics only. To program the instrument please refer to section 6. OPERATING INSTRUCTIONS of the Model 2350-1 Data Logger Manual.

WARNING:

The following procedures require that the instrument be re-initialized. This will clear all logged data and return all programmable settings to their default values. **DETECTOR SETUPS AND LOGGED DATA WILL BE PERMANENTLY LOST.**

1. Turn the Model 2350-1 ON. The display will temporarily go black and then the following two displays will appear for a few seconds each. Then the display that was active when the instrument was turned off will appear.



2. Re-initialize the instrument by entering the cold initialize command (**SSR**).
3. Check the battery voltage on the parameters display (**SVD1**). If it reads below 5.5 volts replace the batteries with fresh alkaline batteries and proceed.

1.3 Threshold and Window Calibration

4. Change the display time base to minutes (**SB1**) display units to counts (**SU7**), and the threshold to 100 (**T100**).
5. Turn the Model 500 Pulser polarity switch to the negative position.
6. Set the pulse amplitude range selector to the 50 mV position, and the LO-HI control to the maximum clockwise position.
7. Set the pulse frequency multiplier to the 1k position and adjust the pulse frequency so that the output is 400k cpm.
8. Connect the Model 2350-1 to the Pulser and verify that the ratemeter reads between 392-408kC/m.
9. Adjust the pulser amplitude control until the meter reads approximately 10 mV. The ratemeter display should read approximately 75% of the count input (**300 kC/m**).

1. INITIAL INSTRUMENT CALIBRATION

10. If the instrument does not read approximately 300 kC/m adjust the gain pot (R10) on the main board until it reads properly.

NOTE: The instrument threshold is now set at 100 = 10 mV

11. Turn the window ON (**WON**), and set it to 100 (**W100**).
12. Adjust the pulse amplitude control until it reads approximately 20 mV. The ratemeter reading should again indicate approximately 75% of the counts (**300 kC/m**).
13. Adjust the threshold setting to 200 (**T200**), 300 (**T300**), and 400 (**T400**) repeating the above procedure each time to verify that the threshold and window settings are linear.
14. Turn the window off (**WOFF**) and reset the threshold to 100 (**T100**).

1.4 High Voltage Calibration

15. Disconnect the instrument from the pulser and connect it to a high impedance voltmeter.
16. Set the high voltage to 1500 volts (**H1500**).
17. If the voltmeter does not read from 1498 - 1502 volts adjust the HV pot (**R52**) on the main board until it reads within tolerance.
18. Adjust the high voltage to 500 volts (**H500**). The voltmeter should read between 490 - 510 volts.
19. Adjust the high voltage to 2000 volts (**H2000**). The voltmeter should read between 1940 - 2060 volts.

1.5 Overload Calibration

21. Disconnect the instrument from the voltmeter, and connect it to the overload simulator resistor box.
22. Change to the Alarm Display (**SVD3**).
23. Turn the overload on (**OON**), and set it to 15.0 μ A (**O150**).
24. Change to the Parameter Display (**SVD1**).
25. Set the high voltage to 1500 volts (**H1500**).
26. Change to the Main Display (**SVD0**). The ratemeter display at the top of the screen should be alternating between the ratemeter display and the word OVERLOAD at approximately 1 second intervals.
27. If this is not occurring adjust the OVERLOAD pot (**R7**) on the main board until the display alternates properly.
28. Adjust the high voltage to 1000 volts (**H1000**). The overload indicator should stop appearing.
29. Adjust the overload to 10.0 μ A (**O100**). The OVERLOAD indicator should again begin to alternate with the ratemeter indicator.
30. Repeat steps 28 and 29 again with a voltage setting of 500 (**H500**) volts and an overload setting of 5 μ A (**O50**).
31. When properly completed turn the overload alarm off (**OOFF**).

MODEL 2350-1 Data Logger

1. INITIAL INSTRUMENT CALIBRATION

1.6 Alarm Setpoints

1. Connect the Model 2350-1 to the pulser.
2. Set the pulse frequency to 600k cpm.
3. Select the Alarm Display (**SVD3**).
4. Set the integrated dose alarm to 100k cpm (**P1E5**).

NOTE: If an alarm is indicated press the ACK/SCRL button to silence the audio.

5. Select the Main Display (**SVD0**) and zero the integrated dose (**SIZ**). The Integrated dose counter should start counting immediately and an alarm should be indicated when it reaches 100k cpm. After the alarm press the ACK/SCRL button to silence the audio.
6. Select the Alarm Display(**SVD3**) and set the scaler alarm to 60,000 counts (**K60000**).
7. Return to the Main Display (**SVD0**), set the scaler count time to 6 seconds (**F6**), and take a scaler count (**C**). When the scaler reaches 60,000 counts an alarm should again be indicated. Press the ACK/SCRL button again to silence the audio.
8. Adjust the pulse frequency of the pulser to less than 600k cpm.
9. Select the Alarm Display (**SVD3**) and set the ratemeter alarm to 600k cpm (**J6E5**).
10. Select the Main Display (**SVD0**) and readjust the pulse frequency to a setting of 600k cpm. The ratemeter alarm should be triggered. Press the ACK/SCRL button to silence the audio.
11. Verify that the ratemeter (**RAT**), scaler (**SCL**), and integrated dose (**DOS**) indicators are all alternating with the alarm (**ALM**) indicators at 1 second intervals.
12. Select the alarm display (**SVD3**)
13. Set all of the alarms to their default high settings (ratemeter, **J1E9**; scaler, **K1000000**; integrated dose, **P1E9**)
14. Select the Main Display (**SVD0**) and reset all of the alarm indicators (**X**).
15. Disconnect the instrument from the pulser and place it into the instrument can.

The electronics have now been calibrated and functionally checked. All remaining calibration procedures involve setting the programmable parameters of the instrument and are covered in the Model 2350-1 Data Logger Instruction Manual.

MODEL 2350-1 Data Logger

2. CALIBRATION ROUTINES

SOURCE PARAMETERS AND TYPICAL VALUES FOR MODEL 2350-1 CALIBRATIONS

To insure that the proper sized sources are used in calculating the calibration constant, and dead time of the Model 2350-1 and detectors the following calculation should be performed.

$$\% \text{ Dead Time Loss} = \text{No. of source counts in 1 second} \times \text{Dead Time (in seconds)}$$

The following list provides some typical LO and HI cal points and dead times for various detectors:

TYPICAL DEAD TIME CALIBRATION, AND OVERLOAD POINTS

DETECTOR	LOW POINT	HIGH POINT	OVERLOAD POINT
44-2	500 μ R/hr	15 mR/hr	100 mR/hr
44-6	20 mR/hr	400 mR/hr	1 R/hr
44-7	5 mR/hr	100 mR/hr	N/A
44-9	6 mR/hr	200 mR/hr	1 R/hr
44-38	20 mR/hr	400 mR/hr	1 R/hr
44-116	6 mR/hr	60 mR/hr	N/A


TYPICAL DEAD TIME AND CALIBRATION CONSTANT SETTINGS

DETECTOR	DEAD TIME	CALIBRATION CONSTANT
44-2	15 - 25 μ sec.	$1.0e^{10} - 0.3e^{10}$
44-6	65 - 135 μ sec.	$7.09e^7 - 0.65e^7$
44-7	250 - 270 μ sec.	$1.39e^8 - 0.1e^8$
44-9	50 - 120 μ sec.	$2.0e^8 - 0.3e^8$
44-38	65 - 135 μ sec.	$7.09e^7 - 0.65e^7$
43-68 Alpha	15 - 25 μ sec.	1
43-68 Beta	18 - 26 μ sec.	1
43-37 Alpha	18 - 20 μ sec.	1
43-37 Beta	19 - 25 μ sec.	1

2. CALIBRATION ROUTINES

2.1 High Voltage Ramp Routine

The high voltage ramp routine is useful when a user needs to plateau a detector to determine the proper operating voltage of the detector. The following display will appear when performing the routine.



HIGH VOLTAGE RAMP ROUTINE DISPLAY

A. INFORMATION PROMPT: Prompts the user for the parameters required for the ramping routine. When the routine is first started the prompt will ask if the user wants to proceed with the current parameters. If answered YES the instrument will then prompt the user to begin the routine by executing the scaler count command. If answered NO then the following information prompts will appear immediately after the one before is answered.

Starting HV ?
Ending HV?
HV Increment?

The message DUMPING HEADER will then appear momentarily followed by the command "Ent C to start".

B. RATEMETER BARGRAPH DISPLAY: A logarithmic display of the ratemeter reading in cps.

C. HV IDENTIFIER: Identifies the current HV setting.

D. STR: Identifies the starting voltage of the HV ramping routine.

E. SCALER COUNT: Shows the scaler reading in progress.

F. SCALER TIMER: Shows the count time for the scaler. When a count is in progress the timer will show the time remaining in the count.

G. END: Identifies the ending HV setting for the ramping routine.

H. INC: Identifies the increment the voltage will be increased by in each step of the routine.

I. OPR CNT: Identifies the count at the beginning knee of the plateau or the recommended operating voltage.

2. CALIBRATION ROUTINES

The following example demonstrates how to run a plateau on a Model 44-2 Gamma Scintillator.

KEYSTROKES

RESPONSE

Note: Before performing the SHR(x) command, certain parameters must be set. The parameters that need to be set include the threshold, window position, and count time.

SHR(x) ENTER Where (x) = the number 0 or 1. 0 = Disables data dump to the serial port. 1 = Enables data dump to the serial port.	The Model 2350 display will show the detector display with the current set of detector parameters. There will be no prompt but you will need to select a set of parameters to use in the routine. If the current parameters are correct press "Y ENTER", otherwise select the set of detector parameters that you want to use and then press "Y ENTER". The instrument will proceed with the routine by displaying the screen shown on the opposite page.
N ENTER	This answers the prompt "PROCEED W/ HV" allowing the user to change the parameters of the routine.
500 ENTER	This set the beginning high voltage at 500 volts.
1000 ENTER	This sets the ending high voltage at 1000 volts.
25 ENTER	This sets the HV increment to 25 volts. (This increment is the amount that the high voltage setting will change between each step of the ramping routine.)
C ENTER	This activates the routine by starting the scaler counting.

"SHR0" will run the plateau routine, determine the operating voltage, but will not dump data to the RS-232 port. "SHR1" will run the plateau routine, determine the operating voltage, and will dump the voltage setting and scaler count at each increment of the plateau routine.


When the routine is complete the following prompt will appear:
"SAVE? HV = XXX"

"Y ENTER" will save the voltage as the operating voltage of the detector and the instrument will return to the screen that was active prior to the execution of the routine. "N ENTER" will return to the screen that was active prior to the execution of the routine without saving the voltage.

2. CALIBRATION ROUTINES

2.2 Calibration Constant and Dead Time Calibration with Background Subtract

This routine will use a two point (low-hi) method to calculate the calibration constant and dead time of the detector. The following display will appear when performing the routine.



CALIBRATION CONSTANT/DEAD TIME ROUTINE DISPLAY

A. INFORMATION PROMPT: Prompts the user for the parameters required for the routine. When the routine is first started the prompt will prompt the user low calibration point. Once entered it will prompt for the following parameters in order.

- HI CAL POINT
- TAKE BACKGROUND
- TAKE LOW SAMPLE
- TAKE HI SAMPLE
- SAVE NUMBERS?

B. RATEMETER BARGRAPH DISPLAY: A logarithmic display of the ratemeter reading in cps.

C. L CAL: Identifies the low calibration point that will be used in the routine.

D. H CAL: Identifies the high calibration point that will be used in the routine.

E. LO SMPL: Shows the reading obtained from the low sample source.

F. HI SMPL: Shows the reading obtained from the high sample source.

NOTE: Line 1 will indicate the ratemeter reading when a sample is being taken, and lines 3 and 4 will show the scaler taking a count. When the routine is complete Line 3 will display the calculated calibration constant, and line 4 will show the calculated dead time constant. These numbers will be saved into the active detector setup if the user answers yes to the prompt to save the numbers.

MODEL 2350-1 Data Logger

2. CALIBRATION ROUTINES

The following example will calculate the cal constant and dead time of a Model 44-2 Gamma Scintillator in R/hr.

KEYSTROKES

RESPONSE

SKB ENTER

The Model 2350-1 will show the calibration constant/dead time routine with background subtract control display with a prompt to set the low cal point.

5e-4 ENTER

This will set the lower calibration point at 500 μ R/hr, then a prompt will appear for the hi cal point.

15e-3 ENTER

This will set the high calibration point at 15 mR/hr, then a prompt will appear to take a background sample.

C ENTER

This will initiate the scaler to take a background reading, then a prompt will appear to take the low sample count. Place a 1 μ Ci ^{137}Cs check source on the end of the detector.

C ENTER

This will initiate the scaler to take a count with the low sample source, then a prompt will appear for the hi sample count. Place a 5 μ Ci ^{137}Cs check source on the end of the detector.

C ENTER

This will initiate the scaler to take the hi sample reading.

NOTE: The SSK command initiates the same routine as stated above except that the background count is not taken or considered when determining the calibration constant and dead time of the system.


ADDITIONAL PARAMETERS

1. The background count should be less than 2% of the low sample count.
2. The low sample count should have at least 15,000 counts.
3. The high sample should be at least 100 times the low sample and produce between 30 - 60% dead time loss.

2. CALIBRATION ROUTINES

2.3 Single Point Calibration Constant Routine with Background Subtract

This routine will use a single point method to calculate the calibration constant of the detector. The following display will appear when performing the routine.



CALIBRATION CONSTANT/DEAD TIME ROUTINE DISPLAY

A. INFORMATION PROMPT: Prompts the user for the parameters required for the routine. When the routine is first started the prompt will prompt the user single calibration point. Once entered it will prompt for the following parameters in order.

TAKE BACKGROUND
 TAKE SAMPLE
 SAVE CAL CNST?

B. RATEMETER BARGRAPH DISPLAY: A logarithmic display of the ratemeter reading in cps.

C. BACKGRND: Identifies the background reading obtained during the routine.

D. SAMPLE: Identifies the sample count taken during the routine.

E. CALPT: Identifies the reference calibration point set for the routine.

NOTE: Line 1 will indicate the ratemeter reading when a sample is being taken, and lines 3 and 4 will show the scaler taking a count. When the routine is complete Line 8 will display the calculated calibration constant. This number will be saved into the active detector setup if the user answers YES to the prompt to save the number.

2. CALIBRATION ROUTINES

The following example will calculate the cal constant of a Model 44-2 Gamma Scintillator in R/hr.

KEYSTROKES

RESPONSE

SSB ENTER

The Model 2350-1 will show the single point calibration constant routine with background subtract control display with a prompt to set the single cal point.

5e-4 ENTER

This will set the single calibration point at 500 μ R/hr, then a prompt will appear to take a background reading.

C ENTER

This will initiate the scaler to take a background count, then a prompt will appear for the sample count. Place the detector in a 500 μ R/hr field.

C ENTER


This will initiate the scaler to take the single sample reading.

NOTE: The SSS command initiates the same routine as stated above except that the background count is not taken or considered when determining the calibration constant.

2. CALIBRATION ROUTINES

2.4 Two Source Dead Time Calibration Routine with Background Subtract

This routine will use a single point method to calculate the calibration constant of the detector. The following display will appear when performing the routine.



CALIBRATION CONSTANT/DEAD TIME ROUTINE DISPLAY

A. INFORMATION PROMPT: Prompts the user for the parameters required for the routine. When the routine is first started the prompt will appear for the background count. Once entered it will prompt for the following parameters in order.

- SAMPLE 1
- SAMPLE 1+2
- SAMPLE 2
- SAVE DEAD TIME?

B. RATEMETER BARGRAPH DISPLAY: A logarithmic display of the ratemeter reading in cps.

C. BACKGRND: Identifies the background count taken during the routine.

D. SAMPLE 1: Identifies the sample count taken with source 1 during the routine.

E. SAMPLE 1 +2: Identifies the sample count taken with sources 1 & 2 during the routine.

F. SAMPLE 2: Identifies the sample count taken with source 2 during the routine.

NOTE: Lines 3 will show a timer and line 4 will show the word COUNTING when a count is in progress. When the routine is complete Line 3 will display the calculated dead time for the detector. This number will be saved into the active detector setup if the user answers YES to the prompt to save the number.

MODEL 2350-1 Data Logger

2. CALIBRATION ROUTINES

The following example will calculate the dead time of a Model 44-2 Gamma Scintillator in R/hr.

KEYSTROKES

RESPONSE

SSD ENTER

The Model 2350-1 will show the dead time calibration routine with background subtract control display with a prompt to take a background reading.

C ENTER

This will initiate a count to get a background sample, then a prompt will appear to take a reading from sample 1.

C ENTER

This will initiate count for sample 1, then a prompt will appear for a count with samples 1 & 2.

C ENTER

This will initiate a count of samples 1 & 2 together, then a prompt will appear for a count with sample 2 only.

C ENTER

This will initiate a count of sample 2 only.

ADDITIONAL PARAMETERS

1. The sources used for this test should be approximately the same size and each should provide a minimum of 25,000 cpm.
2. When combined the sources must provide a between 25 - 60% dead time loss.

3. THEORY OF OPERATION

3.1 AMPLIFIER/POWER SUPPLY BOARD #5371-125

3.1.1 INPUT

Negative going pulses are coupled from the detector through C9 to emitter follower Q1. CR6 is a voltage reference which provides a bias voltage of approximately +3.3 Vdc to Q1 via R45. R9 protects Q1 from input inadvertent shorts/transients. R40 couples the detector to the high voltage supply (HV).

3.1.2 AMPLIFIER

U3A, U3B, U12A, and U12B together are a self-biased amplifier which provides gain in proportion to R12 divided by the series combination of R11 and R10. R10 is adjusted to provide a pulse height in proportion to a given threshold setting at comparator U4 (typically adjusted to trigger with a detector input pulse of 10 millivolts at a threshold setting of 100). Transistor (pins 1, 2, and 6 of U3) provide amplification and pulse polarity inversion. U12A and U12B are a constant current source to pin 6 of U13A. The output self-biases to 2 Vbe — approximately 1.4 Vdc at pin 1 of U13A. This provides just enough bias current through pin 6 to conduct all of the current from the constant current source. Positive pulses from pin 1 are coupled to the Window and Threshold comparators, U4.

3.1.3 WINDOW/THRESHOLD

Threshold (THR) comparator, pins 5, 6, and 7 of U4, provides the lower pulse height discrimination. Pulses are AC coupled from amplifier U3 via C23 to pin 6 of U4. THR reference voltage data from the microprocessor board (μ P) is processed by digital to analog convertor (DAC) U6 and coupled to pin 5 of U4 via opamp U5. R54 and C13 provide additional filtering of the reference voltage. THR reference Vdc is approx. +122 millivolts (mV) with the THR set at 100 on display — 1.222 Vdc at a THR setting of 1000. As the pulse height at pin 6 of U4 increases above the reference

voltage at pin 5, pin 7 — normally high (approx. +5 Vdc) goes low for the pulse duration.

Window (WIN) comparator, pins 1, 2, and 3 of U4 provide the upper pulse height discrimination. Pulses are coupled from amplifier U3 to pin 2 of U4. WIN reference voltage data from the μ P is processed by DAC U6 and coupled to pin 3 of U4. Q2 applies approx. +5 Vdc to the WIN reference to disable the WIN comparator when the WIN OFF signal is applied to the base of Q2. The WIN reference voltage rides on top of the THR reference — i.e., with a THR and WIN setting of 100, the WIN reference at pin 3 of U4 will equal approx. +122 mV referenced above the THR at pin 5 or approx. +244 mV referenced to chassis ground. If the THR is increased to 200 and the WIN still remains at 100 the WIN reference will still equal approx. +122 mV referenced to the threshold but will increase to approx. +366 mV when referenced to ground (THR = 244 mV + 122 mV WIN = 366 mV). As the pulse amplitude increases above the WIN reference voltage, pin 1 of U4 goes low for the pulse duration.

3.1.4 WINDOW/THRESHOLD LOGIC CIRCUIT

Negative pulses from the THR comparator are coupled to univibrator U9. Negative, pulses (approx. 5 volt) are present at univibrator output pin 7 of U9 (PULSE') as long as pins 13 and 3 remain high. When a WIN pulse is present at pins 13 and 3 the Reset function is enabled which disables the PULSE' output locking pin 7 high. Pulses are connected from pin 7 to the μ P on the Central Processor Board for count processing. Pin 7 is tied back to pin 13 via CR1 to provide a time delay (approx. 8 to 10 μ s) for the μ P clock cycle to complete before the next pulse can be recognized by the micro-processor.

3. THEORY OF OPERATION

3.1.5 DIGITAL TO ANALOG CONVERTORS

U6 and U7 are digital to analog convertors (DAC) which convert the digital data from the μ P to analog signals to control the THR, WIN, High Voltage Reference (HV REF), and Overload Reference (OVR REF) variables. Data via BUS0-BUS3 and A0-A2 is loaded into the DAC latches for U6 and U7 by strobing CE3' and CE4' (Chip Select). CE5' (Up Date) is then strobed to transfer the data stored in the latches to the DAC outputs.

3.1.6 HIGH VOLTAGE SUPPLY AND DETECTOR OVERLOAD

Detector High Voltage (HV) is developed by blocking oscillator Q3, T1, C29 and rectified by voltage multiplier CR3-CR5, CR14, CR15, C5, C32, C37, and C38. Q4, CR13, and R41 provide a regulated voltage of approx. 4.4 Vdc (battery voltage must be +4.4 Vdc) to the emitter of Q3. HV increases as current through Q5 increases with maximum output voltage with Q5 saturated.

HV is coupled back through R5, through voltage follower U1, to pin 6 of U2 to complete the regulation loop. Resistors R52 and R6 complete the HV divider network to ground. HV regulation is produced by opamp comparator pins 5, 6, and 7 of U2. During stable operation the voltage at pin 6 will equal pin 5. If HV REF is increased, pin 7 of U2 will increase increasing conduction of Q5 until the voltage at pin 6 equals pin 5 of U2 via HV divider network, R5, R52, and R6. R52 is adjusted to calibrate the HV output to the HV REF signal supplied to pin 5 of U5 by the DAC — HV REF is set at approx. 1500 mV by entering the "H1500" command into the Model 2350-1, the HV CAL is adjusted for 1500 Vdc at the detector connector. C1-C3, C8-C7, C19, and C27 provide additional filtering of the HV output.

Detector Overload is achieved by measuring the voltage differential across R4. As current is increased into the detector, a voltage drop is produced across R4. The HV on either side of R4 is converted to a low voltage by resistor divider networks R5, R52, R6, and R44, R7, and

R43. The voltage differential is coupled to differential opamp, pins 12-14 of U1, via opamp buffers. The differential output is coupled to opamp comparator, pins 8- 10 of U1. OVR REF is provided by the μ P via the DAC — OVERLOAD' is coupled back to the μ P. R7, CURRENT CAL, calibrates the detector current drain to the OVR REF input — i.e., a 10 μ A load is connected to the detector output; the overload is set to "100" by entering the "O100" command (approx. 610mV at OVR REF), R7, CURRENT CAL, is then adjusted until pin 8 or U1 just starts to trigger low.

3.1.7 LOW VOLTAGE SUPPLY

Supply voltages of +5 and -10 Vdc are supplied by U8, T2, and supporting components. U8 is a switching regulator with an on-board comparator providing regulation. An internal voltage reference of 1.0 Vdc maintains the inverting and non-inverting comparators at 1.0V via the feedback loop. Oscillator frequency, set by C36, is approx. 100 kHz at pin 5 of U5. CR11 and C26 provide rectification and filtering for the +5V output. CR12 and C25 provide rectification and filtering for the -10V output.

U11 is a -2.5 Vdc voltage reference for the DAC's — U6 and U7. C8 located on the Backplane board maintains the battery voltage charge during inadvertent battery disconnection (mechanical shock).

3.1.8 LCD VIEWING ANGLE SUPPLY

U2 is configured as a differential opamp which controls the backplane voltage to adjust the viewing angle for the LCD graphics display. Output voltage, pin 1 of U2 is adjustable from 0 to -10 Vdc by varying R50 (VIEWING ANGLE). LCD backplane variance due to temperature is compensated by CR10.

3. THEORY OF OPERATION

3.2 CENTRAL PROCESSOR BOARD #5371-068

3.2.1 MICROPROCESSOR (μ P)

U124, Intel 80C51FA, is the Model 2350-1 central processor. The μ P clock frequency is crystal controlled by Y159 and related components at 6.144 MHz. C172 resets the μ P at power-up to initiate the start of the program routine. The main program routine is stored in EPROM, U122. User parameters and logged data is stored in the 8k X 8 RAM, U123. (replaced with 64k X 8 RAM located on the Memory Expander board)

Address and data information are multiplexed by latch U121 from the μ P to the two memory chips. The low address bits A0-A7 are multiplexed with the data bus, BUS0-BUS7, out of the μ P. U121 latches the address data A0-A7 on lines BUS0-BUS7 during the first part of the external memory cycle when the ALE input is strobed.

The HD175ACCESS shunt configures the μ P parameter access mode and password. When pins 1 and 2 are shunted (default) the access parameters is determined by the settings preset in the EPROM and the password entered in RAM. If shunt is placed between pins 2 and 3 the password is set to "0" and the access level is set to "3".

The ACKNOWLEDGE' input is pulled low when the front panel ACKNL DGE pushbutton switch is depressed silencing the audible alarm(s). The OVERLOAD' signal from the AMP/PS board instructs the μ P to initiate an OVERLOAD condition. WIN OFF' is an output signal from the μ P to the AMP/PS board which disables the Window feature.

The TONE' input comes from the I/O processor which toggles the TONE'/ALARM line low when the external bar code wand registers bar code information. The TONE' and ALARM' (μ P output signal from pin 5 of U124) signals initiate an independent audible tone by changing the CV (control voltage) input on the ICM7556 timer, U118.

The READY, CLEAR SEND, RECEIVE, and XMIT data lines communicate the RS-232 and bar code reader wand information from/to the I/O processor to/from the central processor. The RESET output at pin 13 of U124 is coupled to the

I/O board to reset the I/O processor. Q174 inverts the RESET signal (RESET') providing a reset line for the LCD graphics display. The PULSE' input supplies the pulses from the AMP/PS board to be counted and converted into data by the μ P.

3.2.2 BATTERY VOLTAGE VFC

Un-switched (connected directly to batteries) battery voltage — +BATUNS — is divided by resistor divider network R148 and R149 and connected to voltage to frequency convertor (VFC), U125. VFC OUT is connected to a counter inside the μ P which counts the frequency for a preset time converting to a number displayed on the LCD representing the battery voltage. R130, BAT CAL, adjusts the VFC frequency output to calibrate the LCD BAT reading to the actual battery voltage.

NOTE: Disregard sections 3.2.3 and 3.2.4 for the Model 2350-1

3.2.3 RECORDER DRIVE CIRCUIT

The RCDR output from the μ P is coupled to MOSFET transistor Q126. The μ P uses pulse-width modulation for the recorder signal. The pulse-modulated signal is integrated by R144 and C100. Pins 5, 6, and 7 of U126 "buffers" the integrated voltage from R131. R131 provides calibration of the recorder output voltage. Pins 1, 2, and 3 of U126 are configured as an opamp voltage follower providing recorder drive.

3.2.4 ADDRESS DECODER/VOLTAGE BACKUP

U120 is address decoder which generates the chip enable signals — CEO-CE7 — by monitoring address lines A13-A15. U120 also provides battery backup, +5VBACK, to the Clock chip U119 and the RAM U123. This enables the user to change the instrument batteries without losing RAM memory. +BATUNSR supplies battery voltage to CR112 voltage reference to limit the voltage at 5.1 Vdc. The 5.1 Vdc is coupled to a 0.1

3. THEORY OF OPERATION

farad capacitor through CR113. CR113 prevents C129 from discharging through the +BATUNSR line when the batteries are removed.

3.2.5 CLOCK CHIP

U119 provides the time and date information used in conjunction with the logged data. U119 is crystal controlled by Y158 and supporting components C102 and C103 at 32.786 kHz. The 10 Hz output at pin 13 of U119, provides the timing frequency to the μ P.

3.2.6 AUDIO CIRCUIT

Pulses are coupled from the AMP/PS board to U128 — RDD 104 decade divider chip. AUDIO DIVIDE (front panel) by "1" switch position jumpers the PULSE' input directly to the CLICK' input to the timers. Divide by "10" and "100" switch positions connect the DIV OUT of U128 to the CLICK' input using U128 for pulse division. U118, a dual package 7555 timer, provides audio oscillation. Pulses (CLICK') are coupled to the first timer, pins 1-6 of U118, via C97. This timer is configured as a monostable multivibrator (one-shot) producing a fixed pulse-width for each input pulse at pin 6 of U118. Pulse-width is set by R161 and C95. CR114 limits pulse input to pin 6 of U118 to 5.6 volts; R134 is a pull-up resistor for the AC coupled pulse. The second half of U118 is configured as a stable (free-running) multivibrator which produces the audible tone — frequency is set by R133 and C96.

T157 provides amplification of the audio signal driving the unimorph speaker. Audio volume is varied by adjusting the voltage on top of the primary winding of T157 via the front panel volume control. When an ALARM' is initiated maximum voltage applied by saturating Q154.

3.3 I/O INTERFACE BOARD #5371-063

3.3.1 MICROPROCESSOR (μ P)

U15, Intel 80C51FA, communicates and processes information bi-directionally from/to the RS-232 interface (U13) to/from the central processor via the serial interface lines XMIT, READY, RECEIVE, and CLEAR SEND. U15 also decodes information from the bar code reader (WAND) and sends information to central processor. The μ P clock frequency is crystal controlled by Y159 and related components at 6.144 MHz.

U14, EPROM, stores the program for the I/O μ P. U14 incorporates on-board latches to separate the multiplexed address and data buss lines.

CR26-CR28 provide "clipping" of transients which may be produced at the SERIAL I/O port.

3.3.2 RS-232 INTERFACE

U13, MAX232, is a +5 Vdc powered RS-232 driver/ receiver used to interface the Model 2350-1 to a computer or keypad terminal. U13 incorporates on-board voltage multipliers which generate the +10 and -10 volts required for RS-232 communication.

3.3.3 BACKLIGHT DRIVE

The LCD backlight is illuminated by applying power to +EL POWER when front panel BACKLIGHT switch is ON. With battery voltage applied to +EL POWER, blocking oscillator — Q23, T1 and supporting components begin oscillating. The signal is amplified by T1 to provide a 70VAC signal to drive the LCD backlight. R30, a 10 ohm resistor, limits the power-up surge to the oscillator circuit.

3. THEORY OF OPERATION

3.4 64k MEMORY EXPANDER BOARD #5371-054

NOTE: When the MEMORY EXPANSION board is added, U120 (DS1211) and U123 (CDM6264) on the CENTRAL PROCESSOR board must be removed. The red entries in the CENTRAL PROCESSOR theory of operation should be removed/or inserted for the M2350-1.

U119 is a 128k RAM (only 64k is used) used to store user parameters and logged data. U119 is selected by the central μ P by applying a logic "high" to the RAM input, Pin 30 of U119. The Recorder circuit is disabled on the Central Processor board and the RCDR output on the central μ P is used to enable the Memory Expansion board RAM.

U117 is address decoder which generates the chip enable signals — CE1-CE7, E1 — for the DAC's on the AMP/PS board and the Clock chip on the Central μ P board by monitoring address lines A12-A14 from the central μ P.

U115 provides battery backup, +5VBACK, to the Clock chip (U119) on the Central μ P board and the RAM. This enables the user to change the instrument batteries without losing the RAM memory. U115 also provides a RESET signal which goes low when the +5V supply to pin 2 of U115 drops below a preset threshold (4.5-4.75 volts). This opens U116 analog switches disabling the RAM and CLOCK' signals from the data bus to save the parameter and logging memory.

MODEL 2350-1 Data Logger

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

AMPLIFIER / POWER SUPPLY BOARD 5371-125

Ref. No.	Description	Part No.
CAPACITORS		
C1-C3	0.0056 μ F, 3kV	04-5522
C4	0.1 μ F	04-5521
C5	0.0015 μ F, 3kV	04-5518
C6-C8	0.0056 μ F, 3kV	04-5522
C9	100 pF, 3kV	04-5532
C10	10 μ F	04-5592
C11-C12	100 pF	04-5527
C13	0.01 μ F	04-5523
C14	47 pF	04-5533
C15	100 μ F	04-5576
C16	220 pF	04-5530
C17	100 pF	04-5527
C18-C19	0.0056 μ F, 3kV	04-5522
C20-C21	33 pF	04-5616
C22	0.01 μ F	04-5523
C23	0.1 μ F	04-5521
C24	0.0056 μ F, 3kV	04-5522
C25	22 μ F	04-5579
C26	220 μ F	04-5639
C27	0.0056 μ F, 3kV	04-5522
C28	100 pF, 3kV	04-5532
C29	0.01 μ F	04-5523
C30	100 μ F	04-5576
C31	100 pF	04-5527
C32	0.0015 μ F, 3kV	04-5518
C33	33 pF	04-5616
C34	1 μ F	04-5575
C35-C36	0.001 μ F	04-5519
C37-C38	0.0015 μ F, 3Kv	04-5518
C39	10 μ F	04-5592
C40	0.1 μ F	04-5521
C41	33 pF	04-5616
C42	0.01 μ F	04-5523
DIODES		
CR1-CR2	1N4148	07-6272
CR3-CR5	MR250-2	07-6266
CR6	1N5226	07-6260
CR7	1N5819	07-6306
CR8-CR10	1N4148	07-6272
CR11-CR12	1N5819	07-6306
CR13	1N4148	07-6272

MODEL 2350-1 Data Logger

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

Ref. No.	Description	Part No.
DIODES		
CR14-CR15	MR250-2	07-6266
CR16-CR18	1N4148	07-6272
TRANSISTORS		
Q1	2N3904	05-5755
Q2	MPS6534	05-5763
Q3	MPS-U51	05-5765
Q4	MPS-U01	05-5778
Q5	2N3904	05-5755
RESISTORS		
R1	10kohm, 333MW	12-7748
R2	47 ohm, 333MW	12-7756
R3	56 kohm, 333MW	12-7789
R4	4.7 Mohm	10-7030
R5	1 Gohm	12-7686
R6	470 kohm, 333MW	12-7757
R7	1 Mohm, pot	09-6828
R8	100 kohm, 333MW	12-7747
R9	10 kohm, 333MW	12-7748
R10	100 kohm, pot	09-6823
R11	1 kohm, 333MW	12-7750
R12	470 kohm, 333MW	12-7757
R13-R14	10 kohm, 333MW	12-7748
R15	1 kohm, 333MW	12-7750
R16	4.7 kohm, 333MW	12-7755
R17-R19	10 kohm, 333MW	12-7748
R21	0.1 ohm, 3W	12-7647
R22	80.6 kohm, 333MW	12-7762
R23	1Mohm, 333MW	12-7609
R24	10 Mohm, 333MW	10-7102
R25	10 kohm, 333MW	12-7748
R26	80.6 kohm, 333MW	12-7762
R27	62 kohm, 333MW 1	12-7790
R28	1 Mohm, 333MW	12-7609
R29	33.2 kohm, 333MW	12-7793
R30	47 kohm, 333MW	12-7758
R31	100 ohm, 333MW	12-7746
R32	221 kohm, 333MW	12-7792
R33-R36	10 kohm, 333MW	12-7748
R37	820 ohm, 333MW	12-7791
R38	1 Mohm, 333MW	12-7623

MODEL 2350-1 Data Logger

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

Ref. No.	Description	Part No.
RESISTORS		
R39-R40	1 Mohm	10-7028
R41	100 kohm, 333MW	12-7747
R42	3 kohm, 333MW	12-7768
R43	470 kohm, 333MW	12-7757
R44	1 Gohm	12-7686
R45	100 kohm, 333MW	12-7747
R46	10 kohm, 333MW	12-7748
R47	33.2 kohm, 333MW	12-7793
R48	1 Mohm, 333MW	12-7609
R49	12 kohm, 333MW	12-7787
R51	1.2 kohm, 333MW	10-7116
R52	1 Mohm, pot	09-6828
R53	22 kohm, 333MW	12-7754
R54	10 kohm, 333MW	12-7748
R55	1 kohm, 333MW	12-7750
R56	1 Mohm, 333MW	12-7609
R57	80.6 kohm, 333MW	12-7762
R58	1 Mohm, 333MW	12-7751
R59	10 kohm, 333MW	12-7748
RESISTOR NETWORK		
RN1	5.6 kohm	12-7696
RN20	220 kohm	12-7578
TRANSFORMERS		
T1	M2350-1 HVPS	4275-037
T2	M2350-1 LVPS	4275-089
INTEGRATED CIRCUITS		
U1	LT1079	06-6252
U2	LT1078	06-6251
U3	CMXT3904	05-5888
U4	TLC372	06-6265
U5	LT1078	06-6251
U6-U7	AD7549	06-6253
U8	LM2578	06-6223
U9	CD4098	06-6066
U10	LT1078	06-6251
U11	LM385Z	05-5791
U12	CMXT3906	05-5890
U13	CMXT3904	05-5888

MODEL 2350-1 Data Logger

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

CENTRAL PROCESSOR BOARD 5371-068

Ref. No.	Description	Part No.
CAPACITORS		
C95	0.047 μ F	04-5565
C96	0.001 μ F	14-5519
C97	47 pF	04-5533
C98	0.01 μ F	04-5523
C99-C100	0.1 μ F	04-5521
C101	27 pF	04-5614
C102-C103	22 pF	04-5552
C104	27 pF	04-5614
C109	10 μ F	04-5592
C110	100 μ F	04-5576
C111	1 μ F	04-5575
C129	0.1 F	04-5633
C172	10 μ F	04-5592
C177	0.1 μ F 100V	04-5521
DIODES		
CR112	1N5231	07-6261
CR113-CR115	1N4148	07-6272
CR162	1N4148	07-6272
HEADER		
HD175	3 PIN SIP	n/a
TRANSISTORS		
Q153	2N3904	05-5755
Q154	MPS6534	05-5763
Q155	2N3904	05-5755
Q168,Q174	2N7000	05-5820
RESISTORS		
R130	10 kohm, pot	09-6822
R131	100 kohm, pot	09-6823
R132	2.2 kohm	10-7012
R133	470 kohm	10-7026
R134	220 kohm	10-7066
R135	5.6 kohm	10-7042
R137-R138	10 kohm	10-7016
R139	33 kohm	10-7019
R140	10 kohm	12-7540
R141-R142	100 kohm	12-7557
R143	7.15 kohm	12-7620

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

Ref. No.	Description	Part No.
RESISTORS		
R144	1 Mohm	10-7028
R148	35.7 kohm	12-7640
R149	8.66 kohm	12-7623
R160	220 kohm	10-7066
R161	150 kohm	10-7024
R164	1 Mohm	10-7028
R166-R167	10 kohm	12-7540
R170	1 Mohm	10-7028
R176	220 kohm	10-7066
RESISTOR NETWORKS		
RN151	22 kohm	12-7566
TRANSFORMERS		
T157	M2221, 300-9	4275-074
INTEGRATED CIRCUITS		
U118	ICM7556	06-6244
U119	MM58274	06-6254
U120	CPL-200-T2-T	06-6441
U121	CD74HC573	06-6093
U122	27C512	06-6264
U123	CPL-140-T2-T	06-6439
U124	80C51FA	06-6236
U125	LM331	06-6156
U126	TLC27M71P	06-6248
U128	RDD104	06-6060
CRYSTALS		
Y158	32 kHz XTAL	01-5219
Y159	6.144 MHz XTAL	01-5212

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

64k MEMORY EXPANDER BOARD 5371-054

Ref. No.	Description	Part No.
RESISTORS		
R001	22.1 kohm	12-7843
R114	221 kohm	12-7845
INTEGRATED CIRCUITS		
U115	MAX703ESA	06-6381
U116	CD74HC4066M	06-6323
U117	CD74HC138M	06-6339
U119	CXK581000M	06-6385

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

I/O INTERFACE BOARD 5371-063

Ref. No.	Description	Part No.
CAPACITORS		
C3-C4	27 pF	04-5614
C5	0.047 μ F	04-5565
C6	100 μ F	04-5576
C8-C11	10 μ F	04-5592
C12	4.7 μ F	04-5578
C37-C39	220 μ F	04-5639
RESISTORS		
R16	33 ohm	10-7001
R19	8.2 kohm	10-7015
R20-R21	3.3 kohm	10-7013
R30	10 ohm	10-7046
R40	4.7 ohm	10-7095
RESISTOR NETWORKS		
RN33	10 kohm	12-7777
DIODES		
CR26-CR28	1N5231	07-6261
TRANSISTORS		
Q23	2N3904	05-5755
TRANSFORMERS		
T1	M2350-1 EL	4275-090
INTEGRATED CIRCUITS		
U13	MAX232	06-6188
U36	87C51FA	06-6405
CRYSTALS		
Y25	6.144 MHz MICRO XTAL	01-5212

MODEL 2350-1 Data Logger

4. PARTS LIST, COMPONENT LAYOUTS, AND SCHEMATICS

BACKPLANE BOARD 5371-003

Ref. No.	Description	Part No.
RESISTORS		
R6	1 kohm	10-7009
DIODES		
CR5	1N4148	07-6272
CR7	1N5819	07-6306
CONNECTORS		
J1-J2	125X44-CBNEARS EZA22DRSN	13-8181

See drawings on the following pages.

Model 2350 Amplifier/Power Supply Schematic, Drawing 371 x 125 (3 sheets)

Model 2350 Amplifier/Power Supply Component Layout, Drawing 371 x 127 (2 sheets)

Model 2350-1 Central Processor Board, Drawing 371 x 63

Model 2350-1 Central Processor Board Component Layout, Drawing 371 x 64

Model 2350-1 64K Memory Expander Board, Drawing 371 x 50

Model 2350-1 64K Memory Expander Board Component Layout, Drawing 371 x 51

Model 2350-1 I/O Interface Board, Drawing 371 x 58

Model 2360-1 I/O Interface Board Component Layout, Drawing 371 x 59A


Model 2350-1 Backplane Board, Drawing 371 x 3

Model 2350-1 Backplane Board Component Layout, Drawing 371 x 29


Model 2350-1 Wiring Diagram, Drawing 371 x 72

371125R2P2
 371125R2P2.SchDoc

371125R2P3
 371125R2P3.SchDoc




 LUDLUM MEASUREMENTS, INC.		PO Box 810 501 Oak Street Sweetwater, Texas 79556 U.S.A. 1-800-622-0828		
		Drawn: PAB Design: LL	5/14/2021 5/14/2021	Title: Amplifier/Power Supply Model: 2350 Board#: 5371-125
Approve: <i>[Signature]</i>		Sheet: 1 of 3	Series	Sheet
Print Date: 5/14/2021 1:25:58 PM		Rev: 2	371	125
<small>\\reedont\pcb\Projects\LM\IM 2350\5371-125\Rev2\371125R2P1.SchDoc</small>				




		PO Box 810 501 Oak Street Sweetwater, Texas 79556 U.S.A. +800-822-0828	
		Drawn: PAB Design: LL	5/14/2021 5/14/2021
Board: 5371-125		Sheet: 2 of 3	
Print Date: 5/14/2021 1:25:58 PM	Rev: 2	Series:	Sheet:
U:\temp\Projects\LMIM\2350\5371-125_Rev2\317125R2P1_S48Doc		371	125


A




B



C



D




E

		PO Box 810 501 Oak Street Sweetwater, Texas 79556 U.S.A. 1-800-622-0828	
		Drawn: PAB Design: LL	5/14/2021 5/14/2021
		Board#: 5371-125	
Approve: <i>[Signature]</i> Print Date: 5/14/2021 1:25:59 PM	Sheet: 3 of 3 Rev: 2	Series:	Sheet:
		371	125

\\freedompcb\Projects\LMIM\2350\5371-125\Rev2\371125R2P3.Sch.Doc

Top Layer (Scale 1.3)



 LUDLUM MEASUREMENTS			
Part: 5371-125		Model: 2350	
Desc: Amplifier/Power Supply		Rev:	2
		SHEET: 1 of 4	SERIES SHEET
DWN: PAB	DATE: 5/14/2021	371 127	
APR: <i>PAB</i>	DATE: <i>1/4/2021</i>		

A

B

C

D

Bottom Layer (Scale 1.3)

5371-125

○ ○
REV: 2

8525-040

8525-039

22

1

2

3

4

5

6

1

2

3

4

5

6

A


B

C


D

				LUDLUM MEASUREMENTS			
Part: 5371-125			Model: 2350				
Desc: Amplifier/Power Supply				Rev:	2		
SHEET: 2 of 4			SERIES	SHEET			
DWN: PAB		DATE: 5/14/2021		371	127		
APR: <i>PAB</i>		DATE: <i>5/14/2021</i>					

REVISONS						
EFF	AUTHORITY	ZONE	LTR	DESCRIPTION	DATE	APPROVED




UPDATED	-	LUDLUM MEASUREMENTS INC.				
DR	CKD	26-JUL-99	TITLE: CENTRAL PROCESSOR BOARD			
CHK	RC	7-29-99	BOARD# 5371-868			
USGN	RSS	22-MAR-96	SIZE	MODEL	SHEET	
APPD	RCB	7-29-99	C	2358-1	371	63
NEXT HIGHER ASSY.						
14:08:21		28-Jul-99	SBC71868			SHEET OF 1




LUDLUM MEASUREMENTS INC. SHEETWATER, TX.			
DR	CKB	09/17/96	TITLE: CENTRAL PROCESSOR BOARD
CHK	R.C.	July 28 1999	BOARD: 5371-068
DSGN	RSS	09/17/96	MODEL: 2358-1
APP	251	7-25-99	FILENAME: BS371068
COMPONENT	SOLDER	16:31:18	28-Jul-99
REVISION	SERIES	1.0	371
SHEET			64

- | | |
|---------------|---------------|
| P1 -1 ← NC | P1 -28 ← NC |
| P1 -2 ← A12 | P1 -27 ← NR |
| P1 -3 ← A7 | P1 -26 ← NC |
| P1 -4 ← A6 | P1 -25 ← A8 |
| P1 -5 ← A5 | P1 -24 ← A9 |
| P1 -6 ← A4 | P1 -23 ← A11 |
| P1 -7 ← A3 | P1 -22 ← RD' |
| P1 -8 ← A2 | P1 -21 ← A18 |
| P1 -9 ← A1 | P1 -20 ← NC |
| P1 -10 ← A0 | P1 -19 ← BUS7 |
| P1 -11 ← BUS8 | P1 -18 ← BUS6 |
| P1 -12 ← BUS1 | P1 -17 ← BUS5 |
| P1 -13 ← BUS2 | P1 -16 ← BUS4 |
| P1 -14 ← NC | P1 -15 ← BUS3 |
-
- | | |
|-----------------|-----------------|
| P2 -1 ← VBAT1 | P2 -28 ← +5V |
| P2 -2 ← +SVBACK | P2 -19 ← NC |
| P2 -3 ← NC | P2 -18 ← NC |
| P2 -4 ← NC | P2 -17 ← CE0' |
| P2 -5 ← NC | P2 -16 ← CLOCK' |
| P2 -6 ← NC | P2 -15 ← CE2' |
| P2 -7 ← A15 | P2 -14 ← CE3' |
| P2 -8 ← A14 | P2 -13 ← NC |
| P2 -9 ← A13 | P2 -12 ← CE4' |
| P2 -10 ← GND | P2 -11 ← CE5' |












UPDATED -	/ /	LUDLUM MEASUREMENTS INC			
DR CKD	88/14/95	TITLE: 64K MEMORY EXPANDER BOARD			
CHK	R.C.	BOARD# 5371-854			
DSCH LL	85/12/95	SIZE D	MODEL 2758-1	SERIES 371	SHEET 58
APPD	Jaw	3/7/95			
NEXT HIGHER ASSY.					
14:06:58	26-Feb-98	SBJ7/854		SHEET	UP




LUDLUM MEASUREMENTS INC. SHEETWATER, TX.			
DR	CKB	08/14/95	TITLE: 128K MEMORY EXPANDER BRD
CHK	R.C.	3/4/96	BOARD# S371-054 BS371054
DSGN	LL	05/12/95	MODEL 2350-1 SERIES 371 SHEET S1
APP	JGW	3/7/96	COMP ARTWORK <input type="checkbox"/> SLDR ARTWORK <input type="checkbox"/>
08:38:42	21-Feb-96	COMP OUTLINE <input checked="" type="checkbox"/>	SLDR OUTLINE <input type="checkbox"/>
COMP PASTE <input type="checkbox"/>		COMP MASK <input type="checkbox"/> SLDR PASTE <input type="checkbox"/> SLDR MASK <input type="checkbox"/>	

REVISIONS						
EFF	AUTHORITY	ZONE	LTR	DESCRIPTION	DATE	APPROVED





UPDATED -	/ /	LUDLUM MEASUREMENTS INC.			
DR CHB	11/20/95	TITLE: I/O INTERFACE			
CHK	R.C. 3/4/96	BOARD# 5371-863			
DSCN RSS	11/20/95	SIZE	MODEL	SERIES	SHEET
APPD	JGW 3/7/96	ID	2350-1	371	58
NEXT HIGHER ASSY.					
16:29:49	29-feb-96	SB371063		SHEET	OF 1

371X58




DO NOT PLACE T1 AND R19
PLACE GREEN WIRE JUMPER FROM R19 TO T1


LUDLUM MEASUREMENTS, INC.		PO Box 810 501 Oak Street Sweetwater, TX 79556 U.S.A. 1-800-622-0828	
Title: I/O Interface Board			
Drawn: CKB	11/20/95	Model: 2350	
Design: RSS	11/20/95	Board#: 5371-063	
Approve:		Rev: 1	
PCBA Drawing		SCALE: 1.08	Series Sheet
Print Date: 5/25/2016	1:00:39 PM	Top Overlay	371 59A
W:\Projects\LMM\2350\5371-063\Rev 1\371063R1_Assy.PcbDoc			





UPDATED -	/ /	<input checked="" type="checkbox"/>	LUDLUM MEASUREMENTS INC.
DR LL	9/27/88		TITLE: BACKPLANE
CHK	283 9/10/88		BOARD# 5371-003
DISCN LL	9/27/88		SIZE D
APPD	J6W 9/17/88		MODEL 2350-1
NEXT HIGHER ASSY.			SERIES 371
88:55:18	17-Sep-88	5371-003	SHEET 3 OF 3



<input checked="" type="checkbox"/> LUDLUM MEASUREMENTS INC. SHEETWATER, TX.			
DR	LL	12/17/90	TITLE: BACKPLANE
CHK	R.C.	3496	BOARD# 5J71-003 B5J71003
DSCN	LL	10/20/88	MODEL 2350 SERIES 371 SHEET 29
APP	JGw	3/9/96	COMP ARTWORK <input type="checkbox"/> SLDR ARTWORK <input type="checkbox"/>
08:55:02	21-Feb-96	COMP OUTLINE <input checked="" type="checkbox"/>	SLDR OUTLINE <input type="checkbox"/>
COMP PASTE <input type="checkbox"/>		COMP MASK <input type="checkbox"/>	SLDR PASTE <input type="checkbox"/> SLDR MASK <input type="checkbox"/>



CONNECT PIN 5 TO +BAT UNSWITCHED FOR OPTIONAL #371-067 LASER HAND KIT



DIVIDE 1	DIVIDE 10	DIVIDE 100
1&2	1&2	2&3
4&5	5&6	5&6

UPDATED -		<input checked="" type="checkbox"/> LUDLUM MEASUREMENTS INC.
DR CKB	09/17/96	TITLE: WIRING DIAGRAM
CHK <i>SS</i>	9/17/96	BOARD# 371-071
DSGN RSS	09/17/96	SIZE D
APPD <i>JGw</i>	9/17/96	MODEL 2350-1
NEXT HIGHER ASSY.		SERIES 371
08:53:10	17-Sep-96	4371071
		SHEET 72
		OF 1